**Operating Instructions | Betriebsanleitung | Mode d'emploi |  
Manuale d'istruzioni | Instrucciones de funcionamiento**

# Microsart e.motion

Membrane Filter Dispenser | Membrane Filter-Dispenser |  
Distributeur de membranes filtrantes | Dispensatore di filtri a membrana |  
Dispensador de filtros de membrana

---

<b>English</b>	<b>page 3</b>
<b>Deutsch</b>	<b>Seite 19</b>
<b>Français</b>	<b>page 35</b>
<b>Italiano</b>	<b>pagina 51</b>
<b>Español</b>	<b>página 67</b>

---

# Contents

Intended Use .....	4
Safety Information .....	4
Installation .....	5
Getting Started .....	7
Operation .....	11
Factory Settings .....	13
Troubleshooting .....	14
Maintenance .....	15
Accessories .....	16
Specifications .....	17
Electromagnetic Compatibility .....	18

## Intended Use

The Microsart e.motion dispenser is used to unpack sterile packaged membrane filters conveniently, reliably and quickly and make them easy to remove with tweezers.

The dispenser saves three steps that are otherwise necessary during manual removal, i.e.:

- Opening single packages
- Separating the sensitive membrane
- Removing the membrane from the package

The dispenser can be adjusted so that when tweezers are advanced towards the separator slit it triggers the unpacking of the next membrane. The user can set the sensitivity and the range of the optical motion detector for triggering automatic membrane dispensing.

An optional foot switch is also available.

## Safety Information

The dispenser complies with the European Council Directives as well as international regulations and standards for electrical equipment and electromagnetic compatibility. Improper use or handling, however, can result in damage and/or injury.

- ▶ Read these Operating Instructions thoroughly before using Microsart e.motion

⚠ Do not use Microsart e.motion in zone 0, 1, 20, 21 or 22 hazardous areas|locations.

⚠ If you use electrical equipment in installations and under ambient conditions requiring higher safety standards, you must comply with the provisions specified in the applicable regulations for installation in your country.

- The dispenser may only be used with the following AC adapter: TNG6 Universal, order no. 6971966.
- During cleaning, never allow liquids to enter the dispenser: Only use a slightly moistened cleaning cloth.

▷ If any sign of damage is visible on the equipment or the line cord: Disconnect the equipment from the power supply and make sure that it is protected against further operation.



# Installation

## Equipment Supplied


- The equipment supplied includes the components:
- Dispenser
  - Operating Instructions (this manual)
  - AC adapter with phone plug for connection to the dispenser
  - Line cord with power plug and Euro-compatible IEC plug for connection to the AC adapter

## Connecting to the AC Power


- ▶ Place the dispenser on flat level bench top.


- ▶ Set the power switch (on the back panel of the dispenser) to "0" (middle position).


- ▶ Plug the Euro-compatible IEC plug into the AC adapter.

- ▶ Insert the phone plug into the appropriate socket on the back panel of the dispenser.


- ▶ Plug the power plug into the wall outlet (mains supply).


- ▷ If required, connect the foot switch:  
Connect the plug from the cable on the foot switch into jack on the back panel of the dispenser.

# Getting Started


## Installing the Membrane Filter Strip


- ▶ Set the power switch (on the back panel of the dispenser) to "0" (middle position).


- ▶ Open the membrane cartridge and place it in the dispenser holder so that the cartridge cover can be opened towards the back.


- ▶ Use your left hand to open hinge 1 backwards, and with your right hand feed the membrane filter strips 2 through separator slit 3.

- ▶ Pull membrane filter strips **2** so far forward through the separator slit that **3** and 1/2 packaging units are visible.
- ▶ Flick hinge **1** to the back.
- ▶ Carefully separate the transparent film (top) from the paper carrier strip (bottom).


- ▶ Remove retaining clip **4** from bottom spool **5**.


- ▶ Guide the paper carrier strip around diverter rod **6** and feed roll **7** to bottom spool **5**. Slide retaining clip **5** into the guide slots of the top spool.



- ▶ Wind the paper carrier strip around the bottom spool so that it
  - extends over the two retaining clip guide slots **9**
  - its right edge is up against the guide slot **10** of angular spool.


- ▶ Slide retaining clip **4** into the guide slots on the bottom spool.
- ▶ Keep turning the bottom spool clockwise until the paper carrier strip is under slight pressure.

- ▶ Remove retaining clip **11** from the top spool **12**.

- ▶ Wrap the transparent film around the top spool so that it
  - extends over the two guide slots **13** of the retaining clip
  - its right edge is up against the guide slot of the spool **14**.


- ▶ Slide retaining clip **11** into the guide slots of the top spool.
- ▶ Keep turning the top spool clockwise until the transparent film is under slight pressure.

# Operation


## Selecting an Operating Mode

For dispensing membrane filters, the user has 3 operating modes to select from:

- Transport by pressing the START button
- Transport by pressing the foot switch
- Automatic feed when tweezers are advanced towards the separator slit

- Transport by pressing**
- ▶ Set the power switch (on the back panel of the dispenser) to "I" (top position).
  - > The green LED in the START button lights up continuously.


- Automatic feed**
- ▶ Set the power switch (on the back panel of the dispenser) to "II" (lower position).
  - > The green LED in the START button blinks.

- Transport by pressing the foot switch**
- ▶ If the foot switch is not already connected: connect the foot switch, see also Chapter "Installation".
  - ▶ Set the power switch (on the back panel of the dispenser) to "I" (top position).
  - > The green LED in the START button lights up continuously.


## Unpacking the Next Filter Membrane

### Operating by pressing the start button


- ▶ Press the START button.
- > The membrane filter strip is conveyed one packaging unit further. The next membrane filter is offered for removal.

### Operating using automatic feed

- ▶ Move the tweezers vertically from top to bottom in front of the separator slit.
- > The membrane filter strip is conveyed one packaging unit further. The next membrane filter is offered for removal.

### Operating by pressing the foot switch


- ▶ Operate foot switch.
  - > The membrane filter strip is conveyed one packaging unit further. The next membrane filter is offered for removal.
- 
- ▶ After the last membrane has been used, remove the retaining clips from the spools and pull off the transparent film and paper carrier strip. Slide the retaining clips back on the spools.

# Factory Settings


## Adjusting the Movement Sensor

To automatically dispense membrane filters, the movement sensor for identifying an approaching tweezers can be set at different ranges.

### Low Range Setting

- ▶ Use a screwdriver to turn the potentiometer on the back panel of the dispenser clockwise.
- > Now, bring the forceps closer to the separator slit to trigger the automatic feed mechanism.

### High Range Setting

- ▶ Use a screwdriver to turn the potentiometer on the back panel of the dispenser counterclockwise.
- > The automatic feed mechanism is now triggered when the forceps are at a further distance from the separator slit.


# Troubleshooting

Error	Cause	Solution
2 to 3 filters are feed when button is pressed or motion detected. The red LED in the START button lights up.	The stop sensor in the diverter roller (see figure below) does not identify the rectangular recess in the paper carrier strip.	<b>1)</b> Set power switch to "0" <ul style="list-style-type: none"><li>▶ Detach paper carrier strip from feed roll and spool</li><li>▶ Clean the stop sensor with a dry cloth</li><li>▶ Rewind paper carrier strip</li><li>▶ Set power switch to the desired position</li></ul> <b>2)</b> Sensor test Operate the filter dispenser without any membrane filter strips by pressing the START button. <ul style="list-style-type: none"><li>▶ Sensor OK: The feed roll stops after half a turn.</li><li>▶ Sensor defective: The feed roll does not stop until after approx. 3 sec.</li></ul>

The next filter is not onveyed or only partially when button is pressed or motion detected.  
The red LED in the START button lights up.

Transport of the membrane filter strip is impaired.

▶ Remove feed impairment

The stop sensor is located in the diverter roller which also contains the feed roll.

---

# Maintenance

## **Cleaning the Dispenser**

- △ Make sure that no liquid enters the dispenser
- △ Do not use any aggressive cleaning agents (solvents or similar agents).
- ▶ Disconnect the dispenser from the power supply: remove phone plug from the dispenser.
- ▶ Clean the dispenser using a piece of cloth which has been wet with a mild detergent (soap).
- ▶ After cleaning, wipe down the dispenser with a soft, dry cloth.

# Accessories

Item		Order no.
Foot switch		1ZE---0028
Dispenserfilter:	Contains 3 boxes of 100 membranes	
White with black grid	Ø 47 mm 0.45 µm	11406Z-47----SCM
	Ø 50 mm 0.45 µm	11406Z-50----SCM
High flow	Ø 47 mm 0.45 µm	114H6Z-47----SCM
Grey (black after wetting) with white grid	Ø 47 mm 0.65 µm	13005Z-47----SCM
	Ø 47 mm 0.45 µm	13906Z-47----SCM
	Ø 47 mm 0.8 µm	13004Z-47----SCM
White with green grid	Ø 47 mm 0.45 µm	13906Z-47----SCM
Green with dark-green grid	Ø 47 mm 0.45 µm	13806Z-47----SCM
	Ø 50 mm 0.45 µm	13806Z-50----SCM


# Technical Specifications

Description	Size
Length	204 mm
Height	213 mm
Width	165 mm 250 mm (including removable retaining clip)
Weight	2.9 kg
Operating voltage	110 V or 230 V
Frequency	50 – 60 Hz
Feed time	0.5 s
Operation interrupt after paper feed phase	5 s (prevents inadvertent multiple feed when motion is detected)

**Description of the plug for the foot switch:**

View from the back of the device

---

# Electromagnetic Compatibility

## **CE Mark of Conformity**

The dispenser meets the requirements of the following Directive:  
89/336/EEC "Electromagnetic Compatibility (EMC)"

Applicable European Standards:

Limitation of emissions:

EN 50081-1 Residential, commercial and light industry

EN 50081-2 Industrial environment

Defined immunity to interference:

EN 50082-1 Residential, commercial and light industry

EN 50082-2 Industrial environment

## **Important Note!**

The operator shall be responsible for any modifications to Sartorius Stedim Biotech equipment and for any connections of cables or equipment not supplied by Sartorius Stedim Biotech and must check and, if necessary, correct these modifications and connections. On request, Sartorius Stedim Biotech will provide information on the minimum operating specifications (in accordance with the Standards listed above for defined immunity to interference).

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Phone +49.551.308.0  
Fax +49.551.308.3289  
www.sartorius-stedim.com

Copyright by  
Sartorius Stedim Biotech GmbH,  
Goettingen, Germany.  
All rights reserved. No part of  
this publication may be reprinted  
or translated in any form or by any  
means without the prior  
written permission of Sartorius  
Stedim Biotech GmbH.  
The status of the information,  
specifications and illustrations  
in this manual is indicated by  
the date given below.  
Sartorius Stedim Biotech GmbH  
reserves the right to make changes  
to the technology, features,  
specifications and design of the  
equipment without notice.

Status:  
May 2009,  
Sartorius Stedim Biotech GmbH,  
Goettingen, Germany

Printed in Germany on paper that has  
been bleached without any use of chlorine  
W1A000 · KT  
Publication No.: M--6064-p09054  
Ver. 05 | 2009