




Application lab. report

Nitrided Layer

Preparation to analyze the nitrided layer



QATM-Preparation method

Objective

The aim of the investigation is to determine a suitable preparation method to analyze the nitrided layer for the received samples (Figure 1). The following devices are to use. For cutting the Qcut 400 A, mounting with Qpress 50, grinding and polishing Saphir 560 (we use Saphir 550) and for hardness testing the Qness 60 A+.

Note: The item numbers of the consumables, clamping tools and sample holders could be found in a table at the end of the report.




Figure 1. As-received samples

QATM-Preparation method

Cutting

 Cutting			
Device	Cut-off disc	Anti-corrosion coolant	Clamping tool
Qcut 400 A	Corundum cut-off disc FS-D-350A	QATM Coolant, Standard	- Qtool 80 Vario
Cutting method			
Automatic longitudinal or vertical cutting method (with X or Y-Axis)			
Parameters			
Feed speed	Pulse parameter	Cut-off disc rotational speed	
> 0,3 mm/s	Without pulse	2000 rpm	
Notes			
-			

Hot mounting

 Mounting					
Device	Consumable	Heating time	Temperature	Pressure	Cooling time
Qpress 50	EPO Black	5 Minutes	190° C	300 bar	4 Minutes
Filler or additional consumables	Heating power	Pressure mode	Cooling power		
Bakelit	100 %	Pressure from beginning	100 %		
Notes					
- pre-grind the sample to remove the burr before mounting - this parameter recommendation is for mould ø 40 mm - use two spoon EPO Black and 2 spoon Bakelit					

QATM-Preparation method



Figure 2. Overview of the cutting machine Qcut 400 A



Figure 3. Recommendation for the clamping tool

QATM-Preparation method




Figure 4. Sample before mounting



Figure 5. The mounted sample

QATM-Preparation method

Grinding/Polishing

Device	Samples holder	Pressure mode				
SAPHIR 550	Z5400203	Single				
Step	MEDIUM		 RPM	 N	 min	
 Planar grinding	SiC-Paper, PSA, P320	H ₂ O	200	100 ◀ ▶	15	1:00
 Grinding	SiC-Paper, PSA, P600	H ₂ O	200	100 ◀ ▶	15-20	1:00
 Grinding	SiC-Paper, PSA, P1200	H ₂ O	200	100 ◀ ▶	15-20	1:00
 Polishing	GAMMA	Dia. suspension Wb. Poly, 3 µm	150	100 ▶ ▶	30	5:00
 Polishing	ZETA	Dia. suspension Wb. Poly, 1 µm	150	100 ▶ ▶	30	1:30
 Etching (chem.)	3% alc. NITAL					0:03 – 0:05

Notes

- Pre-dosing for 3 µm und 1 µm: 3 s
Dosing interval and dosing duration for Dia. Suspension 3 µm, 1 µm:
Every 30 s for 1,5 s
- ZETA 1 µm is only to if scratch-free sample surface is needed

QATM-Preparation method




Figure 6. Edge area after GAMMA 3 µm and etching with NITAL – 25:1




Figure 7. Edge area after GAMMA 3 µm and etching with NITAL – 25:1

QATM-Preparation method



Figure 8. Edge area after GAMMA 3 µm - 50:1




Figure 9. Edge area after GAMMA 3 µm and etching with NITAL - 50:1

QATM-Preparation method




Figure 10. Edge area after GAMMA 3 µm and etching with NITAL – 200:1




Figure 11. Edge area after ZETA 1 µm and etching with NITAL – 200:1

QATM- Preparation method

Conclusion

The cutting machine Qcut 400 A can be used very well for the cutting process of the parts. Even if the parts should be bigger it is suitable to use the machine for the cutting process. We recommend the Qtool 80 Vario as a clamping tool. This clamping tool is universal. Big or small samples can be clamped with this clamping tool optimally. We recommend to pre-grind the cutted sample after cutting and before the mounting process.

We recommend to use EPO Black as a mounting resin because this resin has the smallest gap. As a filler-resin please use Bakelit.

Because of the very brittle nitrided surface layer SiC-paper is to use as grinding consumable. To avoid edge rounding we recommend to use the SiC-Paper with self-adhesive (PSA) backside. For easy handling of the self-adhesive SiC-Paper use the Galaxy X-Tap. Grinding till step SiC-paper, P1200 is enough. The polishing step with GAMMA and 3 μm is needed and allows perfect examination of the edge area. The polishing step with ZETA and 1 μm is only to use if absolute scratch-free polishing result is needed.

For etching we recommend 3% alc. based NITAL. Don't do the etching longer the 5 seconds.

The hardness test was made with the hardness tester Qness 60 A+. The test was made with method HV0,1 and HV0,3. Both methods was made on unetched and etched surface. The hardness test report is not included in this report.