Operating Instructions | Mode d'emploi | Instrucciones de manejo

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Cubis®

CUB Models: Semi-microbalances | High-capacity Microbalances

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Contents

1	About These Instructions.....	7
1.1	Validity	7
1.2	Related Documents.....	7
1.3	Target Groups.....	7
1.4	Symbols Used	8
1.4.1	Warnings in Operation Descriptions.....	8
1.4.2	Other Symbols Used.....	8
2	Safety Instructions	9
2.1	Intended Use	9
2.1.1	Modifications to the Device	9
2.1.2	Repairs to the Device.....	9
2.2	Qualifications of Personnel	10
2.3	Functionality of the Device Parts.....	10
2.4	Safety Information on the Device	10
2.5	Electrical Equipment.....	10
2.5.1	Damage to the Device's Electrical Equipment	10
2.5.2	Power Supply Unit and Power Supply Cable.....	10
2.6	Conduct in an Emergency	11
2.7	Accessories and Spare Parts	11
2.8	Risk of Injury During Transporting.....	11
2.9	Glass Breakage.....	11
2.10	Trip Hazard from Connection Cables	11
3	Device Description	12
3.1	Device Overview.....	12
3.2	Draft Shield	13
3.3	Components in the Weighing Compartment	14
3.4	Connections	15
3.4.1	Back of the Device.....	15
3.4.2	Rear Wall of Weighing Compartment	16
3.4.3	Control Unit	16
3.5	Protective Caps and Covers for Connections.....	17
3.6	Conformity-assessed Devices.....	18
3.7	Below-balance Weighing	18
4	Operating Design	19
4.1	Operating Elements in the Main Menu	19
4.2	Operating Elements in Task Management	20
4.3	Status Center	21
4.4	User Management	22
4.4.1	User Profiles	22
4.4.2	User Login	22
4.5	Weighing and Print Profiles	22
4.6	Applications and Tasks	22
4.7	Licensable Device Functions.....	22
4.8	Navigating the Menus	23

5	Installation	25
5.1	Scope of Delivery	25
5.2	Selecting an Installation Site	26
5.3	Unpacking	26
5.4	Attaching or Removing the Control Unit	26
5.5	Placing the Device on its Side and Installing	27
5.6	Inserting the Connection Cable for the Control Unit	27
5.7	Inserting the Connection Cable for Ethernet	28
5.8	Positioning the Weighing Pan and Associated Components	28
5.9	Installing the Draft Shield	29
5.10	Installing the Control Unit	30
5.11	Acclimatizing	30
6	Getting Started	31
6.1	Attaching the Connection Cable for the Control Unit	31
6.2	Connecting the Ethernet Cable	31
6.3	Installing the Power Supply Unit	32
6.4	Connecting the Power Supply	32
6.5	Connecting Accessories	33
6.6	Attaching Protective Caps and Covers	33
7	System Settings	34
7.1	Switching the Device On/Off and Activating Standby Mode	34
7.2	User Login or Logout	34
7.3	Performing System Settings	35
7.4	Using the Help Function	35
7.5	Activating Applications (QAPPs) and Adding to a Task	35
7.5.1	Activating Applications	35
7.5.2	Adding an Application to a Task	36
7.6	Switching Off the isoCAL Function	36
7.7	Activating, Deactivating or Configuring the Ionizer	37
7.8	Configuring Motorized Opening and Closing of the Draft Shield	37
7.9	Adding Weighing and Print Profiles to a Task	38
7.10	Downloading Additional Information	38
7.11	Observing Warm-up Time	38
8	Operation	39
8.1	Manually Opening and Closing the Draft Shield	39
8.2	Opening and Closing the Motorized Draft Shield	39
8.2.1	Opening or Closing by Touching the Touch Sensors	39
8.2.2	Opening or Closing via Proximity Sensors	39
8.2.3	Opening or Closing With the Door Handle	40
8.3	Leveling the Device	40
8.4	Calibration, Adjustment or Linearization	40
8.4.1	Adjusting With the isoCAL Function	41
8.4.2	Internally Calibrating and Adjusting the Device	41
8.5	Preparing Weighings	42
8.6	Weighing	42
8.7	Turning the Ionizer On/Off (Only for Devices with an Ionizer)	43

8.7.1	Setting the Ionizer	43
8.7.2	Starting the Ionization Process	43
8.7.3	Switching Off the Ionizer	43
8.8	Running Application (Example)	43
8.8.1	Executing the "Unit change" Function	43
8.9	Performing Pipette Calibration (Only With Pipette Calibration Kit)	44
8.9.1	Opening and Closing the Closing Slider	44
8.9.2	Filling the Evaporation Trap and Water Reservoir with Liquid	44
8.9.3	Calibrating Pipettes	45
9	Cleaning and Maintenance	47
9.1	Preparing the Device for Cleaning	47
9.2	Cleaning the Device	47
9.3	Recommissioning	47
9.4	Performing a Software Update	48
9.5	Performing a QAPP Center Update	48
10	Malfunctions	50
10.1	Malfunctions on the Operating Display or During Weighing Processes	50
11	Decommissioning	51
11.1	Decommissioning the Device	51
11.2	Dismantling Device Parts	51
11.2.1	Removing the Weighing Pan and Associated Components	51
11.2.2	Dismounting the Draft Shield	52
11.2.3	Dismounting the Connection Cable	53
12	Transport	53
12.1	Transporting the Device	53
13	Storage and Shipping	54
13.1	Storing	54
13.2	Returning the Device and Parts	54
14	Disposal	55
14.1	Disposing of the Device and Parts	55
15	Technical Data	56
15.1	Dimensions and Weights	56
15.2	Ambient Conditions During Storage and Transport	56
15.3	Installation Conditions	56
15.3.1	Installation Site	56
15.3.2	Ambient Conditions at the Installation Site	57
15.4	Acclimatization Before Power Supply	57
15.5	Electrical Data	58
15.5.1	Power Supply	58
15.5.2	Safety of Electrical Equipment	59
15.5.3	Electromagnetic Compatibility	59
15.6	Warm-up Time to Reach Operating Temperature	59
15.7	Interfaces	59

15.7.1	Specifications for the USB-A Interface	59
15.7.2	Specifications for the USB-B Interface	59
15.7.3	Specifications for the USB-C Interface	59
15.8	Model-specific and Licensable Device Functions.....	60
15.9	Ionizer when Using Device in an Isolator with Protective Argon Atmosphere.....	60
15.10	Recommended Calibration Weight.....	60
15.10.1	Semi-microbalance.....	60
15.10.2	High-capacity Microbalance	60
15.11	Conditions for isoCAL Function.....	61
15.12	Data Storage Device	61
15.13	Integrated Clock.....	61
15.14	Backup Battery	62
15.15	Materials	62
15.16	Cleaning Agents and Cleaning Procedures.....	63
15.16.1	Approved Cleaning Agents.....	63
15.16.2	Approved Cleaning Procedures	63
15.17	Metrological Data.....	64
15.17.1	Semi-microbalance.....	64
15.17.2	High-capacity Microbalance	65
16	Accessories and Spare Parts.....	67
16.1	Accessories	67
16.2	Spare Parts.....	69
17	Sartorius Service	70
18	Trademark Information.....	70
19	Conformity.....	70
19.1	EU Declaration of Conformity	70

1 About These Instructions

1.1 Validity

These instructions are part of the device; they must be read in full and kept in a safe place. These instructions apply to the device in the following versions:

Device	Model
Cubis® semi-microbalance	CUB226S5 CUB225S CUB225P CUB125S CUB225S5
Cubis® high-capacity microbalance	CUB116S2 CUB66S CUB66P CUB36S CUB36P CUB26P

1.2 Related Documents

- In addition to these instructions, observe the following documents:
- Instructions for accessories used, e.g., printers or weighing pans
 - Optional: Additional information about cleaning the device (Best Cleaning Practices for Cubis® II Ultra-High Resolution Balances)

1.3 Target Groups

These instructions are addressed to the following target groups. The target groups must possess the specified knowledge.

Target Group	Knowledge and Qualifications
Operator	The operator is familiar with the device and the associated work processes. The operator understands the hazards which may arise when working with the device and knows how to prevent them.*

* If a person in the target group operates the software interface of the device, they are also the “user”.

1.4 Symbols Used

1.4.1 Warnings in Operation Descriptions

WARNING

Denotes a hazard that may result in death or serious injury if it is **not** avoided.

CAUTION

Denotes a hazard that may result in moderate or minor injury if it is **not** avoided.

NOTICE

Denotes a hazard that may result in property damage if it is **not** avoided.

1.4.2 Other Symbols Used

► Required action: Describes activities that must be carried out. The activities in the sequence must be carried out in succession.

▷ Result: Describes the result of the activities carried out.

[] Refers to operating and display elements. Indicates status, warning, and error messages.

 Indicates information for legal metrology for conformity-assessed (verified) devices. Conformity-assessed devices are also referred to as “verified” in these instructions.

Figures in these instructions

Depending on the device configuration, the figures depicting the device and operating display may differ slightly from the supplied device. The variants shown in these instructions are examples.

2 Safety Instructions

2.1 Intended Use

The device is a high-resolution balance, which can be used in laboratories. The device is intended to accurately determine the mass of materials in liquid, paste, powder, or solid form.

Appropriate vessels must be used for loading certain materials.

The device can be operated as follows:

- In individual operation
- Connected to a PC
- Integrated into a network

The device is intended solely for use in accordance with these instructions. Any other use is considered **improper** and can impair the protection supported by the device, e.g., protection against mechanical dangers.

Foreseeable misuse

The following applications are **not** permitted: Operation under any atmosphere other than the normal atmosphere.

Operating conditions for the device

Do **not** use the device in potentially explosive environments. Only use the device indoors.

Do **not** use physical measures to alter the as-delivered state of the device, and only connect approved accessories (see Chapter "16 Accessories and Spare Parts", page 67).

The device may only be used with the equipment and under the operating conditions described in the Technical Data section of these instructions.

2.1.1 Modifications to the Device

If the device is modified: Persons may be put at risk. Device-specific documents and product approvals may lose their validity.

For queries regarding modifications to the device, contact Sartorius.

2.1.2 Repairs to the Device

Device repairs may only be carried out by persons with specialized knowledge of the device. If the device is **not** repaired by a specialist: Persons may be put at risk. Device-specific documents and product approvals may lose their validity.

Sartorius recommends that any repair work, even that carried out after the end of the warranty period, is carried out by Sartorius Service or after consulting with Sartorius Service.

2.2 Qualifications of Personnel

Persons who do not possess adequate knowledge about how to use the device may injure themselves and other persons.

If a particular qualification is required for an activity: The target group is specified. If **no** qualification is specified: The activity can be performed by the "Operator" target group.

2.3 Functionality of the Device Parts

Non-functioning device parts, e.g., as a result of damage or wear, can cause malfunctions. There is a risk of injury to persons.

- ▶ If device parts are **not** functioning: Do **not** use the device.

2.4 Safety Information on the Device

Symbols, e.g., warnings and safety stickers, are safety information for handling the device. Missing or illegible safety information may result in this information **not** being observed. There is a risk of injury to persons.

- ▶ Do **not** conceal, remove, or modify the symbols.
- ▶ Replace the symbols if they become illegible.

2.5 Electrical Equipment

2.5.1 Damage to the Device's Electrical Equipment

Damage to the device's electrical equipment, e.g., damaged insulation, can be life-threatening. Contact with live parts represents a danger to life.

- ▶ If the electrical equipment of the device is defective, disconnect the power supply and contact Sartorius Service.
- ▶ Keep live parts away from moisture. Moisture can cause short circuits.

2.5.2 Power Supply Unit and Power Supply Cable

The use of an **unauthorized** power supply unit or power supply cable may cause life-threatening injuries as a result of electric shocks, for example.

- ▶ Only use the original power supply unit and power supply cable.
- ▶ If the power supply unit or power supply cable needs to be replaced: Contact Sartorius Service. Do **not** repair or modify the power supply unit or power supply cable.

2.6 Conduct in an Emergency

If an emergency occurs, e.g., due to malfunctions of the device or dangerous situations: Persons may be injured. The device must be immediately taken out of operation:

- ▶ Disconnect the device from the power supply.
- ▶ Prevent the device from recommissioning.

2.7 Accessories and Spare Parts

The use of unsuitable accessories and consumables can affect the functionality and operating reliability of the device and has the following consequences:

- Risk of injury to persons
 - Damage, malfunctions, or failure of the device
- ▶ Only use accessories and spare parts that are approved by Sartorius for this device.

2.8 Risk of Injury During Transporting

If the device is **not** transported correctly: The device could fall and cause an injury, e.g., a foot injury. If the device is **not** placed down correctly, e.g., on a laboratory table: Fingers may be crushed.

- ▶ Disconnect the device from all connections at the installation site.
- ▶ Use both hands when transporting the device and setting it down. To do this, reach sideways under the device with both hands in the rear part of the device.
- ▶ Do **not** carry the device on the draft shield or the control unit.

2.9 Glass Breakage

Glass components can break if they fall or are handled incorrectly. Glass fragments can cause cuts.

- ▶ Do **not** use sharp or hard objects on the operating display.
- ▶ **Do not** allow items to fall onto the control unit.
- ▶ In the event of damage to the control unit or draft shield, do **not** use the device. Contact Sartorius Service.

2.10 Trip Hazard from Connection Cables

If the device's connection cables, e.g., the power supply cable, are laid haphazardly: People may trip over the connection cables and injure themselves.

- ▶ Lay all connection cables in such a way that trip hazards are avoided.

3 Device Description

3.1 Device Overview


Fig. 1: Cubis® high-capacity microbalance, with motorized draft shield and 50 mm weighing pan (example)

Pos.	Name	Description
1	Weighing compartment	
2	Draft shield	
3	Cooling element	
4	Weighing module	Contained in the housing.
5	Leveling foot	– For manually leveling the device – Located on the bottom of the device.
6	Proximity sensor	Only high-capacity microbalance: The proximity sensor is intended for contact-free opening or closing of the following components: – Doors of the motorized draft shield – If mounted in the device: Motorized inner draft shield
7	Touch sensor	– Is intended to open or close the following components: – Doors of the motorized draft shield – If mounted in the device: Motorized draft shield – Reacts to touch on the left and right on the top, and on the side panels. – If active: The symbol illuminates at half brightness. – If touched: The symbol illuminates at full brightness.
8	Control unit	Removable.
9	Operating display	Intended to operate the software.
10	Manufacturer's ID label	Contains additional information about the device, e.g., serial number and the metrological and technical data.

3.2 Draft Shield

The draft shield can be opened manually or opening can be motorized, depending on the version of the device.

When opening the draft shield, the top and side doors slide backwards. The top and side doors can be opened together or independently.


Fig. 2: Draft shield, motorized, all doors closed (example)

Pos.	Name	Description
1	Rear panel	Is mounted on delivery.
2	Top door handle	Opens the upper panel.
3	Side door	
4	Bottom door handle	Opens the side door.
5	Front panel	
6	Upper panel	

3.3 Components in the Weighing Compartment


Fig. 3: Components in the weighing compartment and 50 mm weighing pan (example)

Pos.	Name	Description
1	Rear wall of weighing compartment	
2	Guide for glass insert	Available as an accessory.
3	Ionizer nozzle	The ionizer function must be activated in the operating display.
4	Base plate of weighing compartment	Removable, e.g., for cleaning.
5	Pan retainer	Holds the weighing pan. Prevents the weighing pan from turning.
6	Shield plate	Only necessary for 50 mm weighing pan.
7	Weighing pan	50 mm or 90 mm weighing pan (depending on the model)

3.4 Connections

3.4.1 Back of the Device


Fig. 4: Connections on the back of the device, seal for conformity-assessed devices removed

Pos.	Name	Description
1	Slot	For connecting a "Kensington" anti-theft device.
2	USB-C connection	For connecting accessories.
3	USB-A connection	For connecting accessories.
4	USB-B connection	For connection to a PC.
5	On key	Switches the device on.
6	Access switch	Protects the device from changes to the device settings. Is sealed for conformity-assessed devices.
7	Control unit connection	Receives the connection cable for the control unit.
8	Power supply	

3.4.2 Rear Wall of Weighing Compartment


Fig. 5: Connections on the rear wall of the weighing compartment, bezels removed

Pos.	Name	Description
1	Connection for accessories	For connecting accessories, e.g., climate module
2	Connection for accessories	For connecting accessories, e.g., motorized inner draft shield or pipette calibration kit

3.4.3 Control Unit


Fig. 6: Connections to the control unit, some covers removed

Pos.	Name	Description
1	USB-C connection	For connecting accessories.
2	Weighing module connection	For connecting the control unit
3	Ethernet connection	For connection to a network.

3.5 Protective Caps and Covers for Connections


Fig. 7: Protective caps and covers on the rear wall of the weighing compartment, the back of the device, and the control unit

Pos.	Name	Description
1	Cover for rear wall of weighing compartment	Installed in the rear wall of the weighing compartment.
2	Seal for conformity-assessed devices	Screwed onto the device and labeled with a sticker.
3	Protective cap for USB connection	Removable.
4	Cover for below-balance weighing	Screwed on.
5	Protective cap for Ethernet connection	Removable.
6	Cover for battery compartment	Screwed on.
7	Protective cap for USB connection	

3.6 Conformity-assessed Devices

Some settings of conformity-assessed models are protected against operator changes, e.g., “external adjustment”. This measure is intended to ensure the suitability of the devices for use in legal metrology.

3.7 Below-balance Weighing

The device is suitable for below-balance weighing. Samples can be suspended for weighing using below-balance weighing, e.g., samples, which do **not** fit on the weighing pan. Below-balance weighing is possible under the following conditions:

- The device must be set up on a weighing table with recess.
- For below-balance weighing, a below-balance weighing hook must be inserted in the device base. The below-balance weighing hook can be ordered separately as an accessory (see Chapter “16 Accessories and Spare Parts”, page 67).

M

In legal metrology:

- The below-balance weighing equipment may **not** be used.
- The cover of the below-balance weighing may **not** be removed.

4 Operating Design

4.1 Operating Elements in the Main Menu


Fig. 8: Operating elements in the main menu (example)

Pos.	Name	Description
1.	Navigation and function bar	– Enables navigation and searching in menus and lists. – In the “Settings” menu: Displays the name of the menu.
2.	Available tasks	Displays all tasks available for the active user.
3.	Function bar	Displays available submenus and operating functions for the current display and current user.
4.	Task	Starts the described task.

4.2 Operating Elements in Task Management


Fig. 9: Operating elements in task management (example)

Pos.	Name	Description
1.	Navigation and function bar	- Enables navigation and searching in menus and lists. - Enables the addition of tasks. - Opens the QAPP Center. - Displays the name of the menu.
2.	Available tasks	- Displays all available tasks. - Opens a summary of the properties for the displayed task.

4.3 Status Center


Fig. 10: Status center (example)

Pos.	Name	Description
1	Messages	Displays information, warning, and error messages.
2.	Pending transfer	If there are pending transfers: Displays the pending transfers for completed tasks.
3.	Leveling status	Displays the status of the level, starts leveling.
4.	Further categories	Displays further categories that are made visible when swiping to the left, e.g., – Status for the device – Calibration and adjustment report – Audit trail

4.4 User Management

4.4.1 User Profiles

In the factory, four user profiles are created for the device. One role is assigned to each user profile. Each role has rights to operate the device. The rights assigned to each role depend on which device functions the user has to use. User profiles can be adapted.

Additional user profiles, roles, and rights may be activated in the device settings for devices with licensed user management.

4.4.2 User Login

The user must log in to the login display with a user profile. Various setting options and tasks are displayed in the operating display depending on the user profile and role.

4.5 Weighing and Print Profiles

Weighing and print profiles can be created. These profiles can be assigned to a task.

Preset profiles can be used in a task. The preset profiles can be adjusted individually to the application and saved in new profiles for weighing and printing.

4.6 Applications and Tasks

Applications (QAPP applications) are grouped together in QAPP packages. Depending on the model, the device can be supplied with some freely accessible applications. These applications can be used to carry out the most important functions.

Subsequent licensing (package QP99) covers all other applications and can be activated in the QAPP Center for a fee.

An application must be configured as a task before it can be used. Specific settings need to be applied for this, with the help of the assistant function. Tasks are visible to all users who have the required role to perform the task.


4.7 Licensable Device Functions

Some device functions are licensable. In order for these device functions to be usable, they must be activated in the QAPP Center for a fee. The device functions may be activated upon delivery or they can be activated later on in the QAPP Center (licensable device functions see Chapter 15.7, page 59).


4.8 Navigating the Menus

Procedure


- ▶ To open an application or a menu from the main menu:
- ▶ Press the application button or required menu button in the function bar.
- ▶ The application or menu is opened and the name of the open menu is opened in the navigation bar.
- ▶ To return to the main menu from other displays: Press the [Back] or [Menu] button.

- ▶ To return to the next higher menu level: Press the [Back] or [Cancel] button.

- ▶ To scroll through the available tasks or categories in a horizontal menu, e.g., main menu or status center: Swipe to the left or right on the operating display.


- ▶ To scroll through the lists in a vertical menu: Swipe the list downwards or upwards.


- ▶ If a value needs to be selected from a list:
 - ▶ Press the preferred value.
 - ▶ Confirm the selection with [OK].
- ▶ The selected value is saved and the list closes.

- ▶ If elements from a display need to be filtered or a display needs to be browsed:
 - ▶ Press the [Search] or [Filter] button.

- ▶ The keypad is displayed.
- ▶ Type the search value or value to be filtered into the input field (1) using the keypad and confirm with [OK].
- ▶ To close the input field without searching or filtering: Leave the input field empty and press the [OK] button.

- ▶ If language-specific characters need to be entered using the keypad:
 - ▶ Press and hold a letter on the keypad.
 - ▷ If language-specific characters are available for the letter being pressed: The display containing the language-specific characters opens.
 - ▶ Press the required language-specific character.

5 Installation

5.1 Scope of Delivery

Item	Quantity
Device	1
Base plate of weighing compartment	1
Semi-microbalance (models CUB125S CUB225S CUB225P CUB225S5):	
90 mm weighing pan, slotted	1
Semi-microbalance (model CUB226S5) and high-capacity microbalance:	
50 mm weighing pan, slotted	1
Shield plate for 50 mm weighing pan	1
Power supply unit	1
Country-specific power supply cable with test seal	1 – 4*
Control unit with connection cable	1
Draft shield	
Front panel	1
Upper panel	1
Left door	1
Right door	1
In-use dust cover for the control unit	1
Operating Instructions	1 – 2*
Specific certificates for the device	1
* Quantity is country-specific.	

5.2 Selecting an Installation Site

Procedure

- ▶ Ensure that the installation conditions have been met (see Chapter “15.3 Installation Conditions”, page 56).
- ▶ **NOTICE** Risk of damage to the power supply unit from argon! Observe the instructions for using argon (see Chapter “15.3.2 Ambient Conditions at the Installation Site”, page 57).


5.3 Unpacking

The device is packed in multilayer foam packaging. The device components, e.g., the weighing pan, are inserted into the second and third layers of the foam packaging.

The individual layers of the foam packaging must be removed one after the other.

Procedure


- ▶ Open the package.
- ▶ Undo the tension belt.
- ▶ Remove the top layer of foam packaging.
- ▶ Remove the second and third layers of foam packaging.
- ▶ Grip the device by the grooves (1) on the left and right-hand side and lift the device out of the bottom layer of foam packaging.
- ▶ Place the device down on a stable, completely flat surface.
- ▶ Sartorius recommends keeping the original packaging to return the device appropriately, e.g., for repairs.

5.4 Attaching or Removing the Control Unit

Procedure

- ▶ To attach the control unit to the device: Hook the control unit onto the control unit holder (1). The control unit must be positioned entirely on the base.
- ▶ To remove the control unit from the device: Pull the control unit up off of the control unit holder.

5.5 Placing the Device on its Side and Installing

For some installation work, the device has to be placed on its side, e.g., when inserting connection cables.

Material: 1 soft support base, for placing down the device

Requirements

- **No** components are installed in the pan retainer.
- All draft shield doors are removed: Front panel, upper panel, side panels.


⚠ CAUTION

Risk of injury when lifting or transporting!

- ▶ Use both hands when transporting the device and setting it down. To do this, reach sideways under the device with both hands in the rear part of the device.

Procedure


- ▶ Remove the control unit from the device.
- ▶ If the device is to be placed on its side:
 - ▶ Reach sideways under the device with both hands in the rear part of the device.
 - ▶ Turn the device on its side and place it on the soft support base.
- ▶ If the device is to be returned to the upright position:
 - ▶ Reach sideways under the device with both hands in the rear part of the device.
 - ▶ Place the device back onto its base.

5.6 Inserting the Connection Cable for the Control Unit

Procedure

- ▶ Place the device on its side (see Chapter 5.5, page 27).
- ▶ Fully press the connection cable for the control unit into the cable channel (1) along the side of the device. The connection cable for the control unit may **not** come into contact with the side panels as this can cause device malfunctions.
- ▶ Place the device back onto its base.


5.7 Inserting the Connection Cable for Ethernet

An Ethernet connection cable can be attached to the device. The Ethernet connection cable must be inserted into the cable channel on the side of the device.

Material: 1 Ethernet connection cable

Procedure


- ▶ Place the device on its side (see Chapter 5.5, page 27).
- ▶ Fully press the Ethernet connection cable into the cable channel (1) along the side of the device. The Ethernet connection cable may **not** come into contact with the side panels as this can cause device malfunctions.
- ▶ Place the device back onto its base.

5.8 Positioning the Weighing Pan and Associated Components

Procedure


- ▶ Insert the base plate of the weighing compartment (3) into the device.
- ▶ If there is a weighing pan with a shield plate: Place the shield plate (2) onto the pan retainer in the base plate of the weighing compartment.
- ▶ Insert the weighing pan (1) into the pan retainer.
- ▶ If the weighing pan needs to be turned, e.g., with sample holders:
 - ▶ Lift the weighing pan and turn it to the left or right.
 - ▶ The weighing pan snaps into the next position of the pan retainer. This avoids accidental turning of the weighing pan.


5.9 Installing the Draft Shield

Procedure


- ▶ Remove the control unit from the device.
- ▶ Mount the upper panel. To do this, insert the two door panel mounts into the two guide rods on the side of the device and push backwards.

- ▶ Mount the right and left doors. To do this, insert the door panel mounts, one after the other, into the guide rods at the bottom on the side of the device and push backwards.

- ▶ Insert the panel mount for the front panel into the recesses in the front of the device.

5.10 Installing the Control Unit

The control unit can be installed in front of or next to the device.

Procedure

- ▶ Remove the control unit from the device.
- ▶ Place the control unit in the preferred position (dimensions for positioning the control unit see Chapter “15.1 Dimensions and Weights”, page 56). The control unit must be positioned entirely on the base.

5.11 Acclimatizing

When a cold device is brought into a warm environment: The temperature difference can lead to condensation from humidity in the device (moisture formation). Moisture in the device can lead to malfunctions.

Procedure

- ▶ Allow the device to acclimatize at the installation site (acclimatization duration see Chapter 15.4, page 57). The device must be disconnected from the power supply during that time.

6 Getting Started

6.1 Attaching the Connection Cable for the Control Unit

NOTICE


Damage to the device due to incorrect use of the connection cable for the control unit!

If the connection cable for the control unit is connected to other components or to the incorrect connections: The device is damaged.

- ▶ Only use the connection cable for the control unit for the connection between the control unit and the device.

Procedure

- ▶ Insert the connection cable for the control unit into the recess in the rear side of the device (2).
- ▶ Plug the connection cable for the control unit into the "Control unit connection" (1) and screw tightly by hand.


6.2 Connecting the Ethernet Cable

Material: 1 Ethernet connection cable
1 soft surface


Procedure

- ▶ Remove the control unit from the device.
- ▶ Turn over the control unit and place on a soft surface.
- ▶ Open the cable lock (1).
- ▶ Remove the cover for the Ethernet connection from the control unit.
- ▶ Plug the Ethernet cable into the Ethernet connection.
- ▶ Place the Ethernet cable into the cable channel (2) and turn the cable lock (1) using the Ethernet cable.
- ▶ Plug the other end of the Ethernet connection cable into the Ethernet connection at the installation site.


6.3 Installing the Power Supply Unit

Procedure

- ▶ Connect the connection cable from the power supply unit to the “Power supply” connection (1) on the back of the device.

- ▶ Plug the power supply cable into the power supply unit connection.

6.4 Connecting the Power Supply

Requirements

The acclimatization time has been observed and the device has adjusted to room temperature (see Chapter 15.4, page 57).

Procedure

- ▶ Check whether the country-specific power plug matches the power supplies at the installation site.
 - ▶ If required: Contact Sartorius Service.
- ▶ **NOTICE** Damage to the device due to excessive input voltage! Check whether the voltage specifications on the power supply unit match those of the power supply at the installation site.
 - ▶ If the input voltage is too high or too low: Do **not** connect the device to the power supply.
 - ▶ Contact Sartorius Service.
- ▶ Connect the device to the power supply at the installation site. To do this, connect the power plug from the power supply cable into the power socket.
- ▶ The device is switched on and performs the initial functions for device startup.

6.5 Connecting Accessories

Accessories can be connected to the device. If accessories need to be installed in the rear wall of the weighing compartment, e.g., pipette calibration kit: The device must be switched off for the installation process. If accessories are to be connected via a USB connection: The accessories can be connected during operation and installed via the software.

Requirements

The accessories must be suitable for the device (see instructions for the accessories).

Procedure

- ▶ Ensure that the device is disconnected from the power supply.
 - ▶ If required: Disconnect the device from the power supply.
- ▶ If a cover for the rear wall of the weighing compartment needs to be removed from the device: Remove the top (1) and bottom (2) covers for the rear wall of the weighing compartment from the device.
- ▶ Connect the accessories to the appropriate connections on the device (connections for accessory parts see instructions for accessories).
- ▶ Connect the device to the power supply.

6.6 Attaching Protective Caps and Covers

If connections of the device are **not** being used during operation: We recommend using the protective caps and covers provided to seal the connections on the rear wall of the weighing compartment and on the back of the device.

Procedure

- ▶ Check whether all unused connections have been sealed with a protective cap.
 - ▶ If required: Seal the unused connections on the device using the corresponding covers or protective caps.

7 System Settings

7.1 Switching the Device On/Off and Activating Standby Mode


If the device is being connected to the power supply for the first time, or after restoring the factory settings: The device is switched on and the Setup Wizard opens. All steps in the Setup Wizard must be completed.

Requirements

The device is connected to the power supply.

Procedure

- ▶ **NOTICE** Pointed or sharp-edged objects may damage the operating display! Only touch the operating display with your fingertips.
- ▶ If the Setup Wizard is shown: Follow the instructions in the Setup Wizard in the operating display.
- ▶ If the login display is opened: Log into the device using a user profile.
- ▶ To activate standby mode: Press the [On | Off] button.
- ▷ The [Select standby mode] display is opened.
- ▶ Press the required standby mode.
- ▷ The device displays the time.
- ▶ To switch the device off: There are 2 options.
 - ▶ In the [Select standby mode] display, press the [SWITCH OFF DEVICE] button.
 - ▶ Disconnect the device from the power supply.
- ▶ To switch the device back on again from standby mode or after software-controlled switch-off: Press the on key on the back of the device.
- ▶ To switch the device back on again after disconnection from the power supply: Connect the device to the power supply.


7.2 User Login or Logout

The user selection is only displayed when at least one user is created.

Procedure

- ▶ Press the user selection (1).
- ▶ Select a user, e.g., Administrator.
- ▶ Press the [Login] button.
- ▷ If a password is assigned: The input screen for the password opens.
- ▶ Enter the password and confirm.
- ▶ To log out the active user profile from the device: Press the [Logout] button.
- ▶ If required: Log in a new user.

7.3 Performing System Settings

Default settings can be adjusted for the device and the applications in order to align with the ambient conditions and individual operating requirements.

The following settings are necessary to operate the device together with connected components:

- Set up the communication of the connected devices
- Set up additional components

The following settings are recommended to set up the device:

- Set the behavior of the isoCAL function.
- Set the behavior of the motorized draft shield (only for devices with a motorized draft shield).
- If the associated QAPP is activated and the LDAP server is configured: Assign a password.

Procedure

- ▶ Open the main menu.
- ▶ Press the [Settings] button.
- ▶ To adjust settings: Open the preferred submenu.
- ▶ Select the preferred configuration value.
- ▶ Exit the menu.
- ▷ With some settings, the [Booting device] display is displayed in the operating display and the device restarts.

7.4 Using the Help Function

If help texts are available in a menu: The [Help] button is displayed.

Procedure


- ▶ Press the [Help] button.
- ▶ The help texts are displayed.
- ▶ To navigate through the help text: Swipe the text downwards or upwards.

7.5 Activating Applications (QAPPs) and Adding to a Task

7.5.1 Activating Applications

Some applications are activated for the device at the factory. Additional applications may be activated in the QAPP Center. These applications can be tested for 30 days free of charge, and after that require a license.

Procedure

- ▶ Open the task.
- ▶ Press the [QAPP Center] button.
- ▷ An overview of the available QAPP packages is displayed.
- ▶ Select the required QAPP package.
- ▷ A list of the applications contained in the QAPP package is displayed.
- ▶ If the selected QAPP package is to be approved with all the applications it contains:
 - ▶ Press the [License] button.
 - ▷ The input field for the license key is displayed.
 - ▶ If an additional cost is associated with the QAPP package: Enter the license key in the input field and press the [OK] button.
 - ▶ If no additional cost is associated with the QAPP package: Press the [OK] button.
- ▶ If an individual application from the displayed QAPP package is to be activated:
 - ▶ Press the preferred application.
 - ▷ A display opens containing details about the selected application.
 - ▶ Press the [License] button.
 - ▷ The input field for the license key is displayed.
 - ▶ If an additional cost is associated with the application: Enter the license key in the input field and press the [OK] button.
 - ▶ If no additional cost is associated with the application: Press the [OK] button.

7.5.2 Adding an Application to a Task

Applications must be added to a task so that they can run.

Procedure

- ▶ Open the task.
- ▶ Press the [New] button.
- ▷ A list of all activated applications is displayed.
- ▶ To select an application: Press the preferred application.
- ▷ The wizard for creating a new task starts.
- ▶ Follow the wizard's instructions in the operating display.

7.6 Switching Off the isoCAL Function

M

If the isoCAL function is switched off for a conformity-assessed device: The device can only be used in restricted temperature ranges for legal-for-trade applications (see Chapter "15.3.2 Ambient Conditions at the Installation Site", page 57). The isoCAL function **cannot** be switched off for all model versions.

Procedure

- ▶ In the "isoCal execution module" submenu, for the "isoCAL function" parameter, select the "Off" configuration value.

7.7 Activating, Deactivating or Configuring the Ionizer

Requirements

The application for use of the ionizer is activated.

Procedure

- ▶ **NOTICE** Risk of damage to the power supply unit from argon! Observe the instructions for using argon (see Chapter "15.8 Model-specific and Licensable Device Functions", page 60).
- ▶ To deactivate the ionizer: In the menu "Settings" / "Device Settings" / "Ionizer Activation", select the "OFF" configuration value.
- ▶ To activate the ionizer: In the menu "Settings" / "Device Settings" / "Ionizer", select the required intensity and operating time of the ionization process.
- ▶ Exit the settings menu.
- ▷ The [Ionizer] button is displayed in the operating display.

7.8 Configuring Motorized Opening and Closing of the Draft Shield

Activating the touch sensors on the device control unit enables motorized opening and closing of the doors and upper panel of the draft shield.

To activate the touch sensors: Touch the touch sensors or move a hand close enough to the sensor area. The draft shield has a learning capability which enables the following opening parameters to be saved:

- It is possible to control all doors together or individually.
- The opening width of the doors can be adjusted.

Requirements

The application for use of the motorized draft shield is activated.

Procedure

- ▶ Close all draft shield doors.
- ▶ To determine how far a door is opened by using the touch sensors: Manually push the door open to the preferred position.
- ▶ If several doors are to be controlled by the touch sensors simultaneously: Manually push the preferred doors open to the preferred position.
- ▶ Activate the required touch sensor.
- ▷ All open doors are closed.
- ▷ The settings for motorized opening and closing of the draft shield are saved.
- ▷ When the required touch sensor is next touched, the door opens or closes.

7.9 Adding Weighing and Print Profiles to a Task

To be able to use a weighing or print profile: Add a weighing or print profile to a task. Weighing and print profiles can be configured in the Settings menu.

Procedure

- ▶ Open the task.
- ▶ Create or edit a task. To do this, start the wizard to create or edit a task and follow the wizard's instructions in the operating display.

7.10 Downloading Additional Information

As part of the Cubis® IIICUB firmware package on the Sartorius website, additional information is available for the device, e.g., a description of interface protocols, or a set of installation instructions for a website certificate. The information is available as a PDF file, partially in English. The information is stored on the MySartorius website. A Sartorius ID is required to be able to log into MySartorius and call up information. The Sartorius ID can be generated at <https://my.sartorius.com/>.

Procedure

- ▶ Download the "Cubis® CUB Firmware" from the MySartorius website.
- ▶ Find the additional information you require, e.g., the description of interface protocols.

7.11 Observing Warm-up Time

After connecting to the power supply, the warm-up time must be observed. This enables the device to reach its required operating temperature and ensures accurate values during weighing processes.

M

If this relates to a conformity-assessed device: The weight value is marked as **invalid** during the warm-up period.

Procedure

- ▶ Ensure that the warm-up time is observed (see Chapter "15.6 Warm-up Time to Reach Operating Temperature", page 59).

8 Operation

8.1 Manually Opening and Closing the Draft Shield

All doors and the rear panel can be fully or partially opened.

Procedure

- ▶ To open the manual draft shield, e.g., the right-hand door: Take the corresponding door handle and push backwards.
- ▶ To close the manual draft shield, e.g., the right-hand door: Take the corresponding door handle and push it fully forwards.

8.2 Opening and Closing the Motorized Draft Shield

8.2.1 Opening or Closing by Touching the Touch Sensors

Requirements

The motorized opening and closing of the draft shield is configured (see Chapter 7.8, page 37).

Procedure

- ▶ Activate the required touch sensor on the right or left side of the control unit. This enables motorized opening or closing of the motorized draft shield according to the stored setting.

8.2.2 Opening or Closing via Proximity Sensors

The proximity sensor works in "crossover mode":

- Left proximity sensor: Opens and closes the right-hand door and the upper panel
- Right proximity sensor: Opens and closes the left-hand door and the upper panel

The sensitivity of the proximity sensors can be adjusted. Wearing safety gloves may inhibit the use of the proximity sensor.

If a motorized inner draft shield is mounted: The motorized inner draft shield also opens.

Requirements

A license key has been entered to use the motorized draft shield.

Procedure

- Hold your hand in front of the left or right proximity sensor. This enables the door to be fully opened or closed.

8.2.3 Opening or Closing With the Door Handle

Procedure

- To open or close a door: Press on the door handle (1) of a door. The door is automatically fully opened or closed.

8.3 Leveling the Device

Leveling compensates any inclines at the device's installation site. If leveling is necessary: The [Leveling] button is displayed in the weighing display and a message is displayed in the status center.

Procedure

- If the weighing display is opened: Press the [Leveling] button.
- When the status center is displayed: Press the [Level] button.
- ▷ The Leveling Wizard opens.
- Follow the wizard's instructions.

8.4 Calibration, Adjustment or Linearization

Function	Description
Calibration	The device checks how much the displayed value deviated from the specified setpoint.
Adjustment	The device corrects the deviation to the setpoint.
Linearization	The device corrects the deviation from the ideal characteristic curve and the setpoint.

The device needs to be calibrated and adjusted regularly. Various methods can be selected for this:

- Adjusting with the isoCAL function
- Internal calibration and adjustment
- External adjustment
- Internal linearization

Only internal adjustment is described below.

M

External adjustment is **not** possible for conformity-assessed devices in legal metrology.

Procedure

- ▶ Where one of the following conditions occurs, calibrate and adjust the device using the preferred method:
 - Daily, every time the device is switched on
 - After every leveling
 - After changing the ambient conditions (temperature, humidity, or air pressure)
 - After setting the device up at a new installation site

8.4.1 Adjusting With the isoCAL Function

The device can be automatically internally calibrated and adjusted using the isoCAL function.

Requirements

- The isoCAL function is configured in the “Safe weighing” menu, e.g., “On, automatic execution”.
- The conditions for triggering and executing the isoCAL function are met (see Chapter “15.11 Conditions for isoCAL Function”, page 61).

Procedure

- ▶ If the isoCAL function is set to automatic start and the isoCAL function is triggered:
 - ▷ The [isoCAL] button flashes in the operating display.
 - ▶ Wait until the isoCAL function is executed.
 - ▷ In the operating display, a time display counts down from 15 seconds to 0.
 - ▷ If **no** load change or **no** operation takes place on the device before the expiration of the time display: The isoCAL function starts.
- ▶ If the isoCAL function is set to manual start and the isoCAL function is triggered:
 - ▷ The [isoCAL] button flashes in the operating display.
 - ▶ Press the [isoCAL] button.
 - ▷ The isoCAL function starts.
- ▷ If the isoCAL function is complete: The device confirms the completion of the calibration and adjustment process with an acoustic signal, and the calibration report is displayed.
- ▶ To output the calibration report via a connector: Press the [Print Memory] button.
- ▶ To close the calibration report and return to the previous display: Press the [OK] button.

8.4.2 Internally Calibrating and Adjusting the Device

Requirements

The weighing pan is unloaded.

Procedure

- ▶ Open the main menu.
- ▶ Press the “Internal adjustment” task.
- ▷ The internal calibration and adjustment function is executed.
- ▷ If automatic leveling is configured: The device levels itself automatically.
- ▷ If the calibration and adjustment function is complete: The device confirms the completion of the calibration and adjustment process with an acoustic signal, and the calibration report is displayed.
- ▶ To output the calibration report via a connector: Press the [Print Memory] button.
- ▶ To close the calibration report and return to the previous display: Press the [OK] button.

8.5 Preparing Weighings

The device must be prepared before every weighing.

Procedure

- ▶ Level the device.
- ▶ Zero the device. To do this, press the [Zero] button.
- ▶ If the device **cannot** be zeroed: Remove the sample to be weighed and re-zero the device.
- ▶ Adjust the device.

8.6 Weighing

When weighing chemicals, suitable vessels must be used for the material to be weighed. This makes it possible to prevent damage to the device or its accessories.

Requirements

The device is leveled and adjusted.

Procedure

- ▶ Start a task with weighing function.
- ▶ Zero the device. To do this, press the [Zero] button.
- ▶ If below-balance weighing is being carried out: Suspend the sample on the below-balance weighing hook, e.g., with a wire.
- ▶ If using a vessel for the sample:
 - ▶ Place the vessel for the material to be weighed onto the weighing pan.
 - ▶ Press the [Tare] button. This compensates for the weight of the vessel.
 - ▶ Tare the device. To do this, press the [Tare] button.
 - ▶ Place the sample in the vessel or fill the vessel.
- ▶ If **no** vessel is used for the sample: Place the sample on the weighing pan.
- ▶ Once the weight value is displayed in black and the weighing unit is displayed: Read off the measured value.

8.7 Turning the Ionizer On/Off (Only for Devices with an Ionizer)

8.7.1 Setting the Ionizer

Procedure

- ▶ Open the “Settings” / “Device settings” / “Ionizer” menu.
- ▶ Select manual or automatic activation for the “Ionizer Function” parameter.
- ▶ For the “Ionizer intensity” parameter, select the required intensity, e.g., “Weak”.
- ▶ For the “Operating duration” parameter, select the duration of the ionization process, e.g., 60 seconds.

8.7.2 Starting the Ionization Process

Requirements

The model is equipped with an ionizer.

Procedure

- ▶ If the [Ionizer] button is displayed in the weighing display: Press the [Ionizer] button.
- ▷ The ionization process starts.

8.7.3 Switching Off the Ionizer

Procedure

- ▶ Open the “Settings” / “Device settings” / “Ionizer” menu.
- ▶ For the “Ionizer Function” parameter, select the setting value “Off, no function”.

8.8 Running Application (Example)

8.8.1 Executing the “Unit change” Function

The “Unit change” function enables switching between the different units and resolutions defined in the weighing profile of the active task. The units and resolutions can be set at the beginning of the weighing process.

Procedure

- ▶ Start the preferred task.
- ▶ Press the [Unit change] button.
- ▷ All units defined in the weighing profile for the active task are displayed in a list.
- ▷ All resolutions for the weight value defined in the weighing profile for the active task are displayed in a list.

- ▶ Press the preferred unit.
- ▶ To set the resolution for the selected unit: Press the preferred resolution.
- ▶ To confirm the selection and return to the weighing display: Press the [OK] button.
- ▷ The current weight value is displayed in the selected unit and resolution.

8.9 Performing Pipette Calibration (Only With Pipette Calibration Kit)

It is possible to equip the device with the pipette calibration kit accessory and therefore perform pipette calibrations. To do this, the pipette calibration kit must be installed (see the instructions for the accessory).

8.9.1 Opening and Closing the Closing Slider

The closing slider must be opened for pipetting calibration liquid into the pipetting vessel. The closing slider can be opened and closed using gesture control (motion sensor) or the touch sensors. If the closing slider needs to be opened with the touch sensors: The motion sensor is deactivated and is only reactivated when the closing slider has been closed by using the touch sensors.

Procedure

- ▶ To open the closing slider using gesture control: Move hand or pipette tip close to the motion sensor (1).
- ▶ To open the closing slider using the touch sensors: Touch one of the two touch sensors.
- ▷ The closing slider moves to the left and exposes the opening for the pipetting vessel (2).
- ▶ To close the closing slider using gesture control: Move hand or pipette tip away from the motion sensor (1).
- ▶ To close the closing slider using the touch sensors: Touch one of the two touch sensors (3).
- ▷ The closing slider moves to the right and closes the opening for the pipetting vessel.

8.9.2 Filling the Evaporation Trap and Water Reservoir with Liquid

The water reservoir contains the liquid for pipette calibration. Only approved liquids may be used for pipette calibration (approved liquids see the instructions for the accessory). The following information about liquid quantity must be considered when filling the water reservoir:

- If a temperature sensor is used: The temperature sensor must always be covered with liquid during the calibration process.
- The required liquid volume for a complete calibration test series is based on the volume of the pipette to be calibrated and the test conditions according to the ISO 8655 standard.

Material:	– Thermometer
	– Approved liquid, brought to room temperature

Requirements

- The device is leveled and adjusted.
 - If a temperature sensor needs to be used: Temperature sensor is connected (see temperature sensor instructions).
-

NOTICE

Damage to the device due to unacceptable liquid!

If an unsuitable liquid is used for pipette calibration: The device may be damaged. The provisions of the ISO 8655 standard have not been followed.

- ▶ Only use approved liquid (approved liquids see the instructions for the accessory).
-

Procedure

- ▶ Remove the closing slider and the glass cover for the evaporation trap.
- ▶ Pour the approved liquid into the evaporation trap. The bottom must be covered.
- ▶ Reinsert the glass cover and the closing slider.
- ▶ Fill the water reservoir with liquid. Use one of these two options:
 - ▶ Option 1 (**not** applicable with temperature sensor): Remove the water reservoir and fill with approved liquid.
 - ▶ Carefully insert the filled water reservoir into the holder.
 - ▶ Option 2: Pour the approved liquid into the water reservoir.
- ▶ Ensure that the liquids have reached room temperature.
 - ▶ If required: Wait until the liquids in the evaporation trap and the water reservoir have reached room temperature.

8.9.3 Calibrating Pipettes

See the ISO 8655 standard for details about preparing and performing pipette calibration.

Material:	– Exhaust pump
	– Barometer if necessary

Requirements

- The evaporation trap and the water reservoir are filled with enough liquid (see Chapter 8.9.2, page 44.).
- The liquids are at room temperature.

Procedure

- ▶ Prepare pipette calibration in accordance with the requirements in the ISO 8655 standard.
- ▶ Perform pipette calibration in accordance with the requirements in the ISO 8655 standard.
 - ▶ To pipette the absorbed liquid into the pipetting vessel: Open the closing slider (see Chapter 8.9.1, page 44.).

Emptying the pipetting vessel

Procedure

- ▶ Once pipette calibration is completed: Empty the pipetting vessel. To do this, proceed as follows:
 - ▶ Remove the following components one after the other:
 - Closing slider
 - Glass cover
 - Lid
 - ▶ Drain the pipetting vessel with an exhaust pump.

Performing further pipette calibration

Procedure

- ▶ If further pipette calibration needs to be performed:
 - ▶ Re-insert the following device components one after the other:
 - Lid
 - Glass cover
 - Closing slider
 - ▶ Fill the water reservoir with approved liquid (see Chapter 8.9.2, page 44).
 - ▶ Ensure that the requirements have been met.
 - ▶ Prepare and perform pipette calibration again in accordance with the requirements in the ISO 8655 standard.

Finishing pipette calibration

Procedure

- ▶ If **no** further pipette calibration is required:
 - ▶ Empty the evaporation trap with an exhaust pump.
 - ▶ Re-insert the following device components one after the other:
 - Lid
 - Glass cover
 - Closing slider
 - ▶ If a temperature sensor is connected: Remove the temperature sensor and close the opening in the water reservoir.
 - ▶ Empty the water reservoir. Either use an exhaust pump or remove and empty the water reservoir.

9 Cleaning and Maintenance

9.1 Preparing the Device for Cleaning

Procedure

- ▶ If the device is **not** to be cleaned with the Cleaning QAPP: Disconnect the device from the power supply.
- ▶ If an accessory is connected to the device: Disconnect the accessory from the device (see instructions for the accessory).
- ▶ Remove weighing pan and all associated components from the weighing module.
- ▶ If the upper or lower cover for the rear wall of the weighing compartment is removed: Insert the cover for the weighing compartment rear wall into the rear wall of the weighing compartment (see Chapter 6.6, page 33).

9.2 Cleaning the Device

Sartorius recommends cleaning the device at regular intervals, e.g., weekly. There must be **no** foreign substances, e.g., particles, fibers or liquids, in the area of the weighing pan.

To clean the device, you can use Sartorius cleaning utensils or a moistened cloth. The device can be cleaned with the help of the Cleaning QAPP. The Cleaning QAPP contains guided processes for normal and advanced device cleaning. The help texts in the Cleaning QAPP contain notes on valid cleaning agents and cleaning intervals (if cleaning intervals have been set).

Procedure

- ▶  **WARNING** Risk of injury due to electrical current! Protect the power supply unit and the power supply cable from liquids.
- ▶ Only use suitable cleaning agents and cleaning procedures and observe the product information for the cleaning agent used (cleaning agents see Chapter “15.16 Cleaning Agents and Cleaning Procedures”, page 63).
- ▶ If the Cleaning QAPP is to be used for cleaning: Open the Cleaning QAPP for cleaning the device and follow the instructions in the operating display.

9.3 Recommissioning

Procedure

- ▶ Re-insert all components into the device (see Chapter 5, page 25).
- ▶ Connect the required accessories (see Chapter 6.5, page 33).
- ▶ If the device was disconnected from the power supply: Re-connect the device to the power supply (see Chapter 6.4, page 32).

9.4 Performing a Software Update

A software update can be installed (software package) from a USB mass storage device using one of the device's USB connections. An update can also be carried out from a server via other connectors of the device. The following describes the installation of a USB mass storage device, taken from the Sartorius website.

A software update can extend or change the functionality of the device. Sartorius recommends carrying out software updates:

- Prior to starting the software updates, save the device data to a USB mass storage device.
- If a QAPP Center update is to be carried out in addition: Perform the software update for the device first.

Two files are required for the software update: Firmware file with file ending ".upd", and checksum file with file ending ".upd.md5".

Execution of the software update and troubleshooting during this process are described in the help text "Device maintenance".

Requirements

- The device is connected to the power supply.
- The logged-in user has administrator rights.

Procedure

- ▶ Download the software package from the Sartorius website onto the USB mass storage device. To do this, download the "Cubis® III CUB Firmware" file.
- ▶ If this relates to a zip file: Unzip the software package on the USB mass storage device. The files must be saved in the main directory (at root level). The files must **not** be copied into a folder.
- ▶ Insert the USB mass storage device with the software package into the device's USB-A connection.
- ▶ In the "Settings" / "Device maintenance" menu, press the "Update firmware" menu entry.
- ▶ Press the "USB stick" as the connector and select the preferred software version.
- ▷ The software update takes approx. 3 minutes.
- ▷ Once the software update is complete: The software version number is updated in the login display.

9.5 Performing a QAPP Center Update

The QAPP Center package can be installed from a USB mass storage device using one of the device's USB connections. An update can also be carried out from a server via other connectors of the device. The following describes the installation of a USB mass storage device, taken from the Sartorius website.

Sartorius recommends the following when carrying out the QAPP Center update:

- Prior to starting the QAPP Center update, save the device data to a USB mass storage device.
- If a software update is also to be carried out: Perform the software update for the device first.

Two files are required for the QAPP Center update: QAPP Center with file ending “.appcenter”, and checksum file with file ending “qappcenter.upd.md5”.

Execution of the QAPP Center update and troubleshooting during this process are described in the help text “Device maintenance”.

Requirements

- The device is switched on.
- The QAPP Center package is saved on a USB mass storage device or on a server via a connector.

Procedure

- ▶ Download the QAPP Center package from the Sartorius website onto the USB mass storage device. To do this, download the “Cubis® III CUB Firmware” file.
- ▶ If this relates to a zip file: Unzip the QAPP Center package on the USB mass storage device. The files must be saved in the main directory (at root level). The files must **not** be copied into a folder.
- ▶ Insert the USB mass storage device with the QAPP Center package into one of the device’s USB-A connections.
- ▶ In the “Settings” / “Device maintenance” menu, press the “Install QAPP Center” menu entry.
- ▶ Press “USB stick” as the connector.
- ▶ Press the preferred package.
- ▶ Once the QAPP Center update is complete: Confirm successful installation with the [OK] button.
- ▷ Existing tasks remain unchanged after the QAPP Center update. In the existing tasks, the original QAPP versions are used.
- ▶ In order to use the new QAPP version: Create a new task with the new QAPP version. Existing tasks are **not** automatically adapted via the QAPP Center update.

10 Malfunctions

10.1 Malfunctions on the Operating Display or During Weighing Processes

Malfunction	Cause	Solution	Chapter, Page
The operating display is black.	The device is disconnected.	Check whether the power supply unit is connected to the device and the power supply at the installation site.	6.3, 32; 6.4, 32
The operating display is black.	The operating display is not connected.	Check if the connection cable for the operating display is connected to the device.	6.1, 31
The displayed weight value changes constantly.	The device installation site is unstable.	Adjust the parameters in the "Ambient conditions" submenu. Change the installation site.	
	A foreign object is positioned between the weighing pan and the housing.	Remove the foreign object.	
The weight readout displayed by the device is obviously wrong.	The device has not been adjusted.	Adjust the device.	
	The device was not tared before weighing.	Tare the device.	
The operating display is red.	The connection to the device has been interrupted.	Switch the device off and on again.	6.4, 32
An accessory device connected to the USB port is not working.	The connection to the device has been interrupted.	Switch the device off and on again.	7.1, 34
Password forgotten.	A user has forgotten the password.	Contact the administrator to delete or reset the password.	17, 70
	The administrator has forgotten the password.	Contact Sartorius Service to delete or reset the password.	17, 70

11 Decommissioning

11.1 Decommissioning the Device

Procedure


- ▶ Disconnect the device from the power supply. To do this, disconnect the power supply cable from the wall outlet.
- ▶ Disconnect all cables and accessories from the device connections.
- ▶ Place all caps and attachment hoods onto the corresponding connections.
- ▶ If the upper or lower cover for the rear wall of the weighing compartment is removed: Insert the cover for the weighing compartment rear wall into the rear wall of the weighing compartment (see Chapter 6.6, page 33).
- ▶ Clean the device.

11.2 Dismantling Device Parts

11.2.1 Removing the Weighing Pan and Associated Components

Procedure



- ▶ Open the draft shield, e.g., the right-hand door.
- ▶ Remove the weighing pan and the associated components from the device:
 - Weighing pan (1)
 - If there is a weighing pan with a shield plate: The shield plate (2)
 - Base plate of weighing compartment (3)


11.2.2 Dismounting the Draft Shield

Procedure

- ▶ Remove the operating display from the device.
- ▶ Pull off the front panel of the device upwards.

- ▶ Remove the right and left doors. To do this, pull each door forwards to pull it out of the guide rods on the device.


- ▶ Pull the upper panel forwards on both sides and remove it from the guide rods on the device.

11.2.3 Dismounting the Connection Cable

Procedure

- ▶ Place the device on its side (see Chapter 5.5, page 27).
- ▶ Remove the connection cable for the control unit (1) out of the cable channel along the side of the device.
- ▶ Remove the Ethernet connection cable out of the cable channel along the side of the device.
- ▶ Place the device back onto its base.

12 Transport

12.1 Transporting the Device

Requirements

- The device has been decommissioned.
- The control unit is secured to the device.

Procedure

- ▶ **⚠ CAUTION** Risk of injury when lifting or transporting!
 - ▶ Use both hands when transporting the device and setting it down. To do this, reach sideways under the device with both hands in the rear part of the device.
 - ▶ Do **not** carry the device on the draft shield or the control unit.
- ▶ For transporting over longer distances, use a trolley with soft mats. The operating display must be positioned entirely on the base.

13 Storage and Shipping

13.1 Storing

Procedure

- ▶ Decommission the device.
- ▶ Store the device according to the ambient conditions (see Chapter “15.2 Ambient Conditions During Storage and Transport”, page 56).

13.2 Returning the Device and Parts

Defective devices or parts can be returned to Sartorius. Returned devices must be clean and packed in their original packaging.

Transport damage as well as measures for subsequent cleaning and disinfection of the device or parts by Sartorius shall be charged to sender.

Devices contaminated with hazardous materials, e.g., harmful biological or chemical substances, will **not** be accepted for repair or disposal.

Procedure

- ▶ Decommission the device.
- ▶ Contact Sartorius Service for instructions on how to return devices or parts (see www.sartorius.com).
- ▶ Pack the device and its parts in their original packaging for return.

14 Disposal

14.1 Disposing of the Device and Parts

The device and the device accessories must be disposed of properly by disposal facilities.

Two lithium cell batteries, type CR2032, are installed inside the device. Batteries must be disposed of properly by disposal facilities.

Many of the packaging materials are recyclable to promote sustainability and to contribute to reducing global waste volumes.

Procedure

- ▶ Dispose of the device in accordance with local government regulations. Inform the disposal facility that there are two lithium cell batteries, type CR2032, installed inside the device.
- ▶ Dispose of the packaging in accordance with local government regulations.

15 Technical Data

15.1 Dimensions and Weights

	Unit	Value
Dimensions		
Dimensions (L x W x H)	mm	500 x 227 x 302
Dimensions (L x W x H), with control unit removed	mm	365 x 227 x 302
With control unit removed: Max. distance between device and control unit	cm	15
Width of door opening, with fully opened side door	mm	145
Width of door opening, with fully opened upper door	mm	145
Weighing pan dimensions (diameter)		
Semi-microbalance (models CUB125S-3 CUB225S-3 CUB225P-3)	mm	90
Semi-microbalance (model CUB226S-3) and high-capacity microbalance	mm	50
Weight, with draft shield mounted, approx.	kg	14.3

15.2 Ambient Conditions During Storage and Transport

	Unit	Value
Temperature		
During storage and transport	°C	-20 – +60
Storage, dry		

15.3 Installation Conditions

15.3.1 Installation Site

	Unit	Value
Height above sea level, maximum	m	3000
No potentially explosive atmospheres		
Laboratory room, with pollution level according to DIN EN 61010-1		2
Suitable for protection class		
Protection class of the device, according to EN 60529-1		IP30
Protection class of the power supply unit, according to DIN EN 60529		IP30
Access to operation-relevant parts is guaranteed		
Space requirements		

	Unit	Value
Suitable for the dimensions of the device and the associated components		
Setup surface		
Stable, even		
Suitable for the weight of the device and the associated components		
Additional properties		
No heat from heating systems or direct sunlight		
No direct drafts from open windows, AC systems, or doors		
No vibrations		
No "heavy traffic" areas (personnel)		
No electromagnetic fields or electromagnetic radiation, e.g., from radio equipment		
No dry air		
When setting up the device in an isolator with protective argon atmosphere		
Power supply unit is positioned outside the isolator and the argon atmosphere		

15.3.2 Ambient Conditions at the Installation Site

	Unit	Value
Temperature		
In operation	°C	+5 – +40
In operation, with isoCAL function, scope of application as per Directive 2014/31/EU	°C	+10 – +30
In operation, without isoCAL function, scope of application as per Directive 2014/31/EU	°C	+17 – +27
In operation, with conformity-assessed devices, as per specifications on the device ID label		
Relative humidity, in operation		
At temperatures up to 31 °C, maximum	%	80
Then linear decrease from 80 % at 31 °C to 50 % at 40 °C		
Decreasing linearly thereafter, maximum	%	50
Air pressure, minimum	mbar	600

15.4 Acclimatization Before Power Supply

	Unit	Value
Time period between unpacking and connecting to a power supply	h	2

15.5 Electrical Data

15.5.1 Power Supply

	Unit	Value
Power supply only permitted via the power supply cable and power supply unit supplied by Sartorius		
Sartorius power supply unit, type YEPS03-15V0		
Primary (power supply unit)		
AC voltage	V	100 – 240 ± 10 %
Frequency	Hz	50 – 60
Current consumption, maximum	A	1.0
Power consumption, maximum	W	64.5
Secondary (device)		
DC voltage	V	14.25 – 15.75
Current consumption, maximum	A	4.3
Power consumption, typical	W	12
Device fuses		
Quantity		1
Type: Electronic		
Power supply unit fuses		
Quantity		1
Type: Electronic		
Protection class, according to IEC 62368-1		
Device		2
Power supply unit		2
Overvoltage category according to IEC 61010-1		
Device		II
Power supply unit		II

15.5.2 Safety of Electrical Equipment

Safety provisions according to IEC 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements

15.5.3 Electromagnetic Compatibility

Electromagnetic compatibility according to DIN EN/IEC 61326-1: Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 1: General requirements (DIN EN/IEC 61326-1)

Interference immunity: Suitable for use in industrial areas (Table 2 of the standard)

Transient emissions: Class B: Suitable for use in residential areas and areas that are directly connected to a low voltage network that (also) supplies residential buildings

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver
- Consult the dealer or an experienced radio/TV technician for help.

Shielded Cables: Connections between the devices and peripherals must be made using shielded cables in order to maintain compliance with FCC radio frequency emission limits.

Modifications: Any modifications made to this devices that are not approved by Sartorius may void the authority granted to the user by the FCC to operate this equipment.

15.6 Warm-up Time to Reach Operating Temperature

	Unit	Value
Time period between switching the device on and carrying out weighing tasks	h	1

15.7 Interfaces

15.7.1 Specifications for the USB-A Interface

Communication:	USB host (Master)
Connectable devices	Sartorius printers, USB memory sticks, USB barcode scanners, USB keyboards, NFC readers

15.7.2 Specifications for the USB-B Interface

Communication	USB device (Slave)
Type of interface	Virtual serial interface (virtual COM-port, VCP) and "PC-Direct" communication

15.7.3 Specifications for the USB-C Interface

Communication	Downstream-facing port (DFP), USB host (Master)
Communication	RS232 connection with accessory YCC-USB-C-D09M

15.8 Model-specific and Licensable Device Functions

	Semi-microbalance	High-capacity Microbalance
Model-specific device functions		
Proximity sensor	Not available	Available
Licensable device functions		
Use of motorized draft shield	Licensable	Licensable
Use of ionizer	Licensable	Licensable
Use of QAPP-controlled functions for specific tasks on the device, e.g., access management, electronic signature, device cleaning	Licensable	Licensable

15.9 Ionizer when Using Device in an Isolator with Protective Argon Atmosphere

Requirements for the ionizer: Deactivated (in the control unit).

15.10 Recommended Calibration Weight

15.10.1 Semi-microbalance

Models CUB226S5 | CUB225S | CUB225P | CUB125S | CUB125P

		CUB226S5	CUB225S	CUB225P	CUB125S	CUB225S5
	Unit	Value	Value	Value	Value	Value
Resolution	mg	0.005	0.01	0.01 0.1	0.01	0.05
Supported load	g	220	220	120 220	120	220
Recommended accuracy class		E2	E2	E2	E2	E2
External test weight	g	200	200	200	100	200

15.10.2 High-capacity Microbalance

Models CUB116S2 | CUB66S | CUB66P

		CUB116S2	CUB66S	CUB66P
	Unit	Value	Value	Value
Resolution	mg	0.002	0.001	0.01 0.001
Supported load	g	111	61	61 12
Recommended accuracy class		E2	E2	E2
External test weight	g	100	50	50

Models CUB36S | CUB36P

		CUB36S	CUB36P	CUB26P
	Unit	Value	Value	Value
Resolution	mg	0.001	0.01 0.001	0.01 0.001
Supported load	g	32	32 10	32 10
Recommended accuracy class		E2	E2	E2
External test weight	g	20	20	20

15.11 Conditions for isoCAL Function

	Unit	Value
Possible conditions for triggering the isoCAL function		
In the event of a temperature change	K	1.5
After a time interval	h	12
After successful leveling		
Only conformity-assessed devices: After interrupting the power supply		
Required conditions for executing the isoCAL function		
Operating display is in weighing mode (not in menu)		
Alphanumeric inputs are not active		
Min. time period without entry on the device	min	2
Min. time period with unchanged loading of the weighing pan	min	2
Max. loading of weighing pan, high-capacity	%	2

15.12 Data Storage Device

	Value
Maximum number of entities	500000

15.13 Integrated Clock

	Unit	Value
Maximum deviation per month (RTC)	s	30

15.14 Backup Battery

	Unit	Value
Lithium battery, type CR2032		
Service life at room temperature, minimum	years	10

15.15 Materials

Housing

Stainless steel 1.4401 | 1.4404, Aluminum

Plastic PBT | PA

Float glass Optiwhite

Control unit

Plastic PA12

Operating display: Plastic PA12 | float glass

Weighing pan

Top section: Titanium

Bottom section: Stainless steel, zinc oxide, brass

15.16 Cleaning Agents and Cleaning Procedures

15.16.1 Approved Cleaning Agents

Device Components	Cleaning Agents and Concentration					
	Ethanol, 70 %	Isopropanol, 70 %	Citric acid, 10 %	Diluted hy- drogen per- oxide, 3.5 %	Sodium hy- droxide, 32 %	Ecolab Klarcide™ Sporicidal Chlorine Thiosulfate
Draft shield						
Top door	x	x	xx	xx	-	xx
Left and right doors	x	x	xx	xx	-	x
Front panel	x	x	x	x	-	x
Components in the weighing compartment						
Weighing pan	x	x	x	x	xx	x
Shield plate	x	x	x	x	xx	x
Base plate of weighing compartment (removable)	x	x	x	x	xx	x
Rear wall of weighing compartment	xx	x	x	x	x	x
Weighing compartment support (holds the base plate of the weighing compartment)	x	x	x	x	x	x
Control unit	x	x	x	x	x	x
Back of the device						
Plastic surfaces	x	xx	x	x	x	x
Cooling element	x	xx	x	x	x	x
x Suitable xx Suitable, may affect visual appearance of the device, does not impact mechanical stability - Not suitable						

15.16.2 Approved Cleaning Procedures

Wiping the device surfaces with a lightly moistened cleaning cloth		
Spraying the device surfaces with cleaning agents, contact time	min	5 – 10
Drying the surfaces with a non -abrasive cloth.		

15.17 Metrological Data

15.17.1 Semi-microbalance

Models CUB226S5 | CUB225S | CUB225P | CUB125S | CUB125P

		CUB226S5	CUB225S	CUB225P	CUB125S	CUB225S5
	Unit	Value	Value	Value	Value	Value
Scale interval (d)	mg	0.005	0.01	0.01 0.1	0.01	0.05
Maximum capacity (Max)	g	220	220	120 220	120	220
Repeatability at up to 5 % load						
Standard deviation of the load values, tolerance	mg	0.01	0.015	0.015	0.015	0.04
Standard deviation of the load values, typical value	mg	0.004	0.0065	0.0065	0.0065	0.02
Repeatability at approx. maximum capacity						
Standard deviation of the load values, tolerance	mg	0.025	0.025	0.04	0.025	0.06
Standard deviation of the load values, typical value	mg	0.015	0.015	0.02	0.015	0.025
Linearity deviation						
Tolerance	mg	0.07	0.07	0.1	0.07	0.1
Typical value	mg	0.03	0.03	0.03	0.03	0.05
Deviation when load is off-center, positions according to OIML R76						
Test weight	g	100	100	100	50	100
Tolerance	mg	0.12	0.15	0.2	0.12	0.2
Typical value	mg	0.04	0.05	0.06	0.04	0.1
Sensitivity drift of +10 °C to +30 °C	ppm/K	1	1	1	1	1
Tare maximum capacity: Less than 100 % of maximum capacity						
Accuracy class, according to Directive 2014/31/EU		I	I	I	I	I
Verification scale interval (e), according to Directive 2014/31/EU	mg	1	1	1	1	1
Minimum load (Min), according to Directive 2014/31/EU	mg	1	1	1	1	1
Minimum initial weight according to USP (United States Pharmacopeia), Chap. 41						
Optimum minimum initial weight	mg	4.1	8.2	8.2	8.2	41
Typical stabilization time	s	1.5	1.5	1.5	1.5	1.5
Typical measurement time	s	6	4	4	4	4

15.17.2 High-capacity Microbalance

Models CUB116S2 | CUB66S | CUB66P

		CUB116S2	CUB66S	CUB66P
	Unit	Value	Value	Value
Scale interval (d)	mg	0.002	0.001	0.001 0.01
Maximum capacity (Max)	g	111	61	12 61
Repeatability at up to 5 % load				
Standard deviation of the load values, tolerance	mg	0.004	0.0015	0.002
Standard deviation of the load values, typical value	mg	0.0025	0.0007	0.0007
Repeatability at approx. maximum capacity				
Standard deviation of the load values, tolerance	mg	0.01	0.004	0.01
Standard deviation of the load values, typical value	mg	0.005	0.0025	0.006
Linearity deviation				
Tolerance	mg	0.03	0.02	0.02
Typical value	mg	0.02	0.005	0.008
Deviation when load is off-center, positions according to OIML R76				
Test weight	g	50	20	20
Tolerance	mg	0.03	0.02	0.03
Typical value	mg	0.02	0.01	0.012
Sensitivity drift of +10 °C to +30 °C	ppm/K	1	1	1
Tare maximum capacity: Less than 100 % of maximum capacity				
Accuracy class, according to Directive 2014/31/EU		I	I	I
Verification scale interval (e), according to Directive 2014/31/EU	mg	1	1	1
Minimum load (Min), according to Directive 2014/31/EU	mg	0.2	0.1	0.1
Minimum initial weight according to USP (United States Pharmacopeia), Chap. 41				
Optimum minimum initial weight	mg	1.64	0.82	0.82
Typical stabilization time	s	3.5	3.5	3.5 2.5
Typical measurement time	s	8	10	10 6

Models CUB36S | CUB36P

		CUB36S	CUB36P	CUB26P
	Unit	Value	Value	Value
Scale interval (d)	mg	0.001	0.001 0.01	0.001 0.01
Maximum capacity (Max)	g	32	10.1 32	5.1 22
Repeatability at up to 5 % load				
Standard deviation of the load values, tolerance	mg	0.0015	0.002	0.002
Standard deviation of the load values, typical value	mg	0.0007	0.0007	0.0007
Repeatability at approx. maximum capacity				
Standard deviation of the load values, tolerance	mg	0.0025	0.007	0.007
Standard deviation of the load values, typical value	mg	0.0018	0.005	0.005
Linearity deviation				
Tolerance	mg	0.012	0.015	0.015
Typical value	mg	0.005	0.006	0.006
Deviation when load is off-center, positions according to OIML R76				
Test weight	g	10	10	10
Tolerance	mg	0.015	0.02	0.02
Typical value	mg	0.006	0.008	0.01
Sensitivity drift of +10 °C to +30 °C	ppm/K	1	1	1
Tare maximum capacity: Less than 100 % of maximum capacity				
Accuracy class, according to Directive 2014/31/EU		I	I	I
Verification scale interval (e), according to Directive 2014/31/EU	mg	1	1	1
Minimum load (Min), according to Directive 2014/31/EU	mg	0.1	0.1	0.1
Minimum initial weight according to USP (United States Pharmacopeia), Chap. 41				
Optimum minimum initial weight	mg	0.82	0.82	0.82
Typical stabilization time	s	3.5	3.5 2.5	3.5 2.5
Typical measurement time	s	10	10 6	10 6

16 Accessories and Spare Parts

16.1 Accessories

This table contains a selection of the accessories that can be ordered. For information on other products, contact Sartorius.

Item	Quantity	Order Number
Thermal transfer direct thermal printer for GLP GMP printouts on continuous paper and tags	1	YDP30 YDP30-NET
Wireless nano USB adapter for company network or independent Wi-Fi network, e.g., operation with a Sartorius network printer YDP30-NET (only for Europe)	1	YWLAN01MS
Printer paper and thermal transfer ribbon, 90 m	1	69Y03285
Self-adhesive paper and thermal transfer ribbon, 90 m	1	69Y03286
Thermal paper, 24 m	5	69Y03287
Self-adhesive thermal paper, 24 m	5	69Y03288
Tag roll, 58 x 30 mm	1000	69Y03092
Tag roll, 58 x 76 mm	500	69Y03093
Tag roll, 58 x 100 mm	350	69Y03094
Wireless nano router, e.g., for Sartorius network printer YDP30-NET for operation in an independent Wi-Fi network (only for Europe)	1	YWLAN02MS
Connection cable for operating display, length 3 m	1	YCC01-CUB-3
Dust cover	1	YDCC2CUB
Connection cable with RS232 adapter, USB-C to RS232, 9-pin	1	YCC-USB-C-D09M
Ethernet extension, 1 m	1	YCC-RJ45-CAT7
Motion sensor with USB connection cable	1	YHS02USB
Foot switch for draft shield, tare, print	1	YFS02
Inner draft shield		
Motorized	1	YDS125A
Manual	1	YDS125U
Glass insert to reduce the height of the weighing compartment	1	YDSHR
Climate module (connection to rear wall of weighing compartment)		
For determining temperature, humidity and air pressure	1	YCM20MC
For determining temperature, humidity and air pressure, with DAkkS calibration certificate	1	YCM20MC-DAKKS
Climate module tower (USB connection), for determining temperature, humidity and air pressure	1	YCM20MC-TOWER
Glass insert to reduce the height of the weighing compartment	1	YDSHR
Rear wall with syringe holder for infusion pump tests	1	YIP-3

Item	Quantity	Order Number
Sample holder		
Configurable for sample bulbs of up to 50 mL	1	YSH02-3
For Safe-Lock® 2 mL sample tubes	1	YSH14-3
For Safe-Lock® 5 mL sample tubes	1	YSH18-3
For bottles up to 40 mL	1	YSH22-3
For coronary stents	1	YSH12-3
For syringes	1	YSH46-3
Filter weighing pan, including the shield plate for rear wall of weighing compartment		
For filters up to 75 mm	1	YSH35-3
For filters up to 150 mm	1	YSH30-3
Weighing scoop holder	1	YSH26-3
Weighing pan		
90 mm weighing pan, slotted	1	YWP10-3
50 mm weighing pan, slotted, with shield plate for 50 mm weighing pan	1	YWP09-3
Dust cover	1	6960MS09
Below-balance weighing hook for below-balance weighing	1	69MS0307
Density determination kit for solids and liquids, models: 116S2, 226S-3, 125S, 225P, 225S, 225S5	1	YDK03MC
Pipette calibration kit for semi-microbalances and high-capacity microbalances	1	YCP07MC
Cubis® pipette calibration kit	1	YCP01MPS
Weighing table		
Weighing table made from synthetic stone, with vibration dampening	1	YWT03
Weighing table made from wood with synthetic stone	1	YWT09
Wall console	1	YWT04
Weighing scoop made from chrome-nickel steel, (90 mm x 32 mm x 8 mm)	1	641214
Cleaning kit	1	YCK01MC
RFID reader with USB-A connection	1	YRFID01
USB QR barcode scanner	1	YBR05
Ionizing blower	1	YIB01-ODR
Ionizer including universal adapter	1	YIB03-C
Sartorius ionizer Stat-Pen	1	YSTPP01
Software: QAPP licenses		
Package with applications for ionizer	1	QAPP1001

Item	Quantity	Order Number
Package with applications for motorized draft shield	1	QAPP1002
Package with applications for hardware	1	QP10
Package with all applications	1	QP99

16.2 Spare Parts

This table contains a selection of the spare parts that can be ordered. For information on other products, contact Sartorius.

Item	Quantity	Order Number
Left door for draft shield (external draft shield)	1	YCCDSL
Right door for draft shield (external draft shield)	1	YCCDSR
Upper door for draft shield (external draft shield)	1	YCCDSU
Front panel for draft shield (external draft shield)	1	YCCDSF
Power supply	1	69MS0320

17 Sartorius Service

Sartorius Service is available for queries regarding the device. Visit the Sartorius website (www.sartorius.com) for information about the service addresses, services provided, or to contact a local representative.

For inquiries about the system or when contacting Sartorius Service in the event of a malfunction, keep the device information close at hand, e.g., serial number, hardware, firmware, configuration, to pass it on to Sartorius Service. To do this, refer to the information on the manufacturer's ID label and in the "General device Information" menu.

18 Trademark Information

Ecolab Klercide™ is a registered trademark of Ecolab Europe GmbH.

19 Conformity

19.1 EU Declaration of Conformity

The attached Declaration of Conformity hereby confirms compliance of the device with the directives cited.

The Declaration of Conformity supplied with the balance is for conformity-assessed (verified) balances for use in the EEA. Keep it in a safe place.

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Device type

***Semi micro and high-capacity Micro Balance | Power Supply | Hand sensor | Climatic
module | Automatic motorized inner draft shield***

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*in the form as delivered fulfils all the relevant provisions of the following European Directives
and meets the applicable requirements of the harmonized European Standards including any
amendments valid at the time this declaration was signed listed below:*

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The person authorised to compile the technical file: Sartorius Lab Instruments GmbH & Co. KG
Electronics & Product Compliance
37070 Goettingen, Germany

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Phone: +49 551 308 0
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote all genders.

Copyright notice:

These instructions, including all components, are protected by copyright.

Any use beyond the limits of the copyright law is not permitted without our approval.

This applies in particular to reprinting, translation and editing irrespective of the type of media used.

Last updated:

04 | 2026

© 2026

Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Str. 20
37079 Goettingen, Germany

SP | Publication No.: WMC6035-g260402