



# EECPC CALIBRATION: ORDER THE RIGHT EQUIPMENT

---

APPLICATION NOTE EECPC-005 (US)

Calibration of the Engine Exhaust Condensation Particle Counter (EECPC), TSI model 3790A, requires following the procedures described in EURO 5 Regulation 83. Using these systems these components adhere to the guidelines of ISO 27891.

The instruments needed to construct the calibration bench are similar to standard particle instruments from TSI® Incorporated. See Figure 1 for a schematic of the entire calibration process.

---

## Aerosol Generation

- 3074B Pressurized Air Supply. This optional item is used to assure that the pressurized air used in the calibration process is dry, and free of oil and particles.
- 3480 Electro spray Aerosol Generator (EAG)
- 1050001 Dilution Bridge
- 3082 Electrostatic Classifier
- 3088 Soft X-ray Neutralizer (strongly recommended; in general, it increases the concentration of particles available for calibration)
- 3085A Nano DMA

---

## Flow Balance

- 1602051 Filter
- 4148 Flow Calibrator (includes small filter shown in Figure 1) (permits calculation of relative humidity of air presented to aerosol instrumentation, when ambient humidity is known)

---

## EECPC Calibration

- 3708 Flow Splitter
- 3790A Engine Exhaust CPC (CPC under Test)
- 3750 Condensation Particle Counter (Reference CPC)
- 3068B Aerosol Electrometer (Reference Aerosol Electrometer)

---

## Tubing

- 3001788 Conductive Tubing (fits ¼" barb)
- 3001789 Conductive Tubing (fits ⅜" barb)

Details for generating emery oil aerosol for the calibration process are provided in a separate application note.





**Figure 1: Schematic of the EECPC Calibration Setup**



TSI Incorporated – Visit our website [www.tsi.com](http://www.tsi.com) for more information.

USA	Tel: +1 800 680 1220	India	Tel: +91 80 67877200
UK	Tel: +44 149 4 459200	China	Tel: +86 10 8219 7688
France	Tel: +33 1 41 19 21 99	Singapore	Tel: +65 6595 6388
Germany	Tel: +49 241 523030		