

Findenser™ Mini NS24 Cone, NS24 Socket

Technical datasheet



Weight	670 g
Length × Diameter	275 mm × 72 mm
Joints	B14, B19, B24
Inner or wetted parts	Borosilicate glass
Finned jacket	Anodised aluminium
Plug & Seal	Acetal & FKM
Thermofluid	Water
Max solvent bpt	155 °C
Oven drying	Max Temp. 60 °C
Dishwasher	Max Temp. 50 °C
Autoclaving	Do Not Autoclave
Ambient operating temp.	0 °C to 24 °C