



Findenser™ NS24 Cone, NS24 Socket

Technical datasheet

Weight	1200 g
Length × Diameter	400 mm × 72 mm
Joints	B19, B24, B29
Inner or wetted parts	Borosilicate glass
Finned jacket	Anodised aluminium
Plug & Seal	Acetal & FKM
Thermofluid	Water
Max solvent bpt	155 °C
Oven drying	Max Temp. 60 °C
Dishwasher	Max Temp. 50 °C
Autoclaving	Do Not Autoclave
Ambient operating temp.	0 °C to 24 °C