

Hei-FLOW Core 120 Multi

Technical datasheet



Dimensions (W × H × D)	166 × 225 × 256 mm, w/o pump head
Weight	7.6 kg (w/o pump head)
Drive	EC motor
Speed range	10-120
Display	-
Protection class (EN 60529)	IP 54
Acoustic pressure	< 55 (dB(A)) (in accordance with IEC 61010)
Rated voltage	110-240 V, 50/60 Hz
Protection class	1 (EN 61140)
Overheat protection	yes
Control accuracy	± 0,5 %
Speed control	digital
Interfaces	-
Power consumption	100 W
Drive power	100 W
Operation mode	continuous 24/7
Direction of rotation	left / right
Turns indicator	analogue
Volumetric flow single channel (ml/min)	0.85-861
Volumetric flow multi channel (ml/min)	0.005-364
Precision	±5 %

Volume dosing (ml)	-
Interval dosing (ml)	-
Soft start	-
Electronic brake	-
Foot switch port	-