

Hei-VOLUME Distimatic Pro Industrial with adaptor for A-series glassware

Technical datasheet



Dimensions (W × H × D)	186 × 429 × 521 mm, including control panel
Weight	approx. 17 kg
Protection class (EN 60529)	IP42
Rated voltage	100 – 240 V AC; 50/60 Hz
Connection	L+N+PE
Protection class	I
Overvoltage category	II
Degree of pollution	2
Power input	1,500 W max
Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 m asl
Fuse	15 A