

Reax 2 incl. universal adaptor

Technical datasheet



Dimensions (W × H × D)	501 × 235 × 180 mm
Weight	approx. 5 kg
Maximum permissible load	1 kg
Drive	Shaded pole motor
Speed range	20 – 200 rpm; Indicator: analogue
Protection class (EN 60529)	IP21
Acoustic pressure	< 45 (dB(A)) (in accordance with IEC 61010)
Rated voltage	230-240 V (50 Hz) or 115 V (60 Hz)
Protection class	I
Overvoltage category	n.a.
Degree of pollution	n.a.
Power input	33 W max.
Operating temperature	n.a.
Installation altitude	n.a.
Overheat protection	self reset function
Speed control	electronic
Usable area	230 × 160 mm
Movement	Overhead rotation
Orbit	n.a.
Fuse	n.a.
Permissible mains voltage fluctuations	n.a.

Bottle dimensions

n.a.

Interfaces

n.a.