

Unimax 1010

Technical datasheet



Dimensions (W × H × D)	320 × 125 × 375 mm
Weight	approx. 8 kg
Maximum permissible load	5 kg
Drive	EC motor
Speed range	30 – 500 rpm
Protection class (EN 60529)	IP30
Acoustic pressure	< 50 dB (A), in accordance with IEC 61010
Rated voltage	100 – 240 V, 50/60 Hz
Protection class	I
Overvoltage category	n.a.
Degree of pollution	n.a.
Power input	50 W
Operating temperature	5 °C – 31 °C up to 80 % rel. humidity 32 °C – 40 °C up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 m asl
Overheat protection	n.a.
Speed control	n.a.
Usable area	290 × 258 mm
Movement	orbital
Orbit	10 mm
Fuse	M 1,25 A (230 V) or T 2,0 A (115 V)

Permissible mains voltage fluctuations	n.a.
---	------

Bottle dimensions	n.a.
--------------------------	------

Interfaces	1 × RS232 (Sub-D, 9 pole)
-------------------	---------------------------