

SYSTEM LAB FURNITURE

HIGH-QUALITY STANDARD CABINETS

The concept of our system lab furniture is founded on 4 corner stones:

QUALITY

Criteria for our furniture are based on the exacting requirements of industry and laboratories. The combination of high-quality materials and carefully thought-out construction guarantees stability and lasting quality.

FLEXIBILITY

The QATM-system provides a customized and flexible platform for personal workspaces, entire buildings, and everything in between.

FUNCTIONALITY

Perfect fit with QATM equipment and their installation. Archiving solutions and covered cable guides ensure a well-organized and efficient workplace.

DESIGN

Clean lines and high-quality materials create a fascinating, timeless and homogeneous style.

HIGHLIGHTS



Ideal basis for placing QATM benchtop machines



Sturdy, corrosion-free construction, powder-coated aluminum housing




Fronts and worktops made of water-resistant HPL solid material



Ideal for integration of recirculation cooling system or settling tank



Soft closing doors and drawers



PREMIUM QUALITY
MADE IN GERMANY



Easy transportation using a pallet truck



Integrated cable and media duct



Installation area for media connections on the back side

part of **VERDER**
scientific

AVAILABLE STANDARD CABINETS

CABINET DIMENSIONS	W X H X D	SUITABLE FOR
2-part framework	880 x 900 x 800 mm	for all QATM benchtop machines, except Qpol twin wheel units, except cut-off machine Qcut 350 A
3-part framework	1320 x 900 x 800 mm	for all QATM benchtop machines, for Qcut 350 A use LAB CABINET #9



LAB CABINET #1

- 2-part framework
- 2 doors,
- 1 utility port (large brush)

Order No. M5810001



LAB CABINET #2

- 2-part framework
- 2 doors, 2 drawers,
- 1 utility port (large brush)

Order No. M5810002



LAB CABINET #3

- 2-part framework
- 1 door, 3 drawers
- 1 utility port (large brush)

Order No. M5810003



LAB CABINET #4

- 2-part framework
- 1 door, 5 drawers,
- 1 utility port (large brush)

Order No. M5810004



LAB CABINET #5

- 2-part framework
- 2 doors,
- 1 utility port (large brush)
- 1 sink, ceramic

Order No. M5810005



LAB CABINET #6

- 3-part framework
- 3 doors,
- 2 utility ports (large brush, cover)

Order No. M5810006



LAB CABINET #7

- 3-part framework
- 3 doors, 3 drawers
- 2 utility ports (large brush, cover)

Order No. M5810007



LAB CABINET #8

- 3-part framework
- 2 doors, 3 drawers
- 2 utility ports (large brush, cover)

Order No. M5810008



LAB CABINET #9

- 3-part framework
- 2 doors, 4 drawers
- 2 utility ports (large brush, cover)

reinforced version

Order No. M5810009



LAB CABINET #10

- 3-part framework
- 2 doors, 5 drawers
- 2 utility ports (large brush, cover)

Order No. M5810010



LAB CABINET #11

- 3-part framework
- 3 doors,
- 2 utility ports (large brush, cover)

1 sink, ceramic

Order No. M5810011




LAB CABINET #12

- corner framework,
- 1 corner door,
- 2 utility ports (large brush, cover)

Order No. M5810012

AVAILABLE ACCESSORIES


With the right accessories, you have the opportunity to further optimize the functionality and efficiency of your cabinets.



FOR GRINDING PAPER

- disc storage for grinding paper, with 6 compartments

Order No. A5810352



FOR GRINDING DISCS

- disc storage for grinding discs, with 13 compartments

Order No. A5810351



DRAWER INSERT SETS

- insert set 1 set 2

Order No. Z5810001 Z5810002

- insert set 3 set 4

Order No. Z5810003 Z5810004

