



pushing boundaries

Linseis Messgeräte GmbH | Vielitzerstr. 43 | D-95100 Selb

Specifications UDSC L64

	UDSC L64
Temperature range	From -60 °C to 160°C
Ramp speed	0,001 to 10 °C/min
Control accuracy	100µ°C
Sample volume	5 to 100µL
Sensitivity	550µV/mW
Unique features	Direct T sample measurement Joule effect calibration
Equipment dimensions	L*w*h = 700*650*500 mm

UST-Id.Nr.DE 132 951 806 | TELEFON +49 9287 880-0 | FAX +49 9287 70488 | E-MAIL: info@linseis.de

Bankkonto: Sparkasse Hochfranken IBAN DE 89 7805 0000 0760 0044 65 BIC BYLADEM1HOF

Bankkonto: Commerzbank IBAN DE 39 7004 0041 0228 8256 00 BIC COBADEFFXXX

Geschäftsleitung: Dipl. Phys. Claus Linseis • Registergericht: HRB HOF ABT. B NR: 299